

Datasheet: 9100-0004

Description:	SHEEP ANTI HUMAN TRANSFERRIN
Specificity:	TRANSFERRIN
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Human
-----------------------	-------

Product Form	Purified IgG - liquid
---------------------	-----------------------

Antiserum Preparation Antisera to human transferrin were raised by repeated immunisations of sheep with highly purified antigen. Purified IgG prepared by affinity chromatography on Protein G

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
---------------------------------	--

Approx. Protein Concentrations	IgG concentration 5.0 mg/ml
---------------------------------------	-----------------------------

Immunogen	Pure human transferrin prepared from pooled normal human serum.
------------------	---

External Database Links

UniProt:
[P02787](#) [Related reagents](#)

Entrez Gene:[7018](#) TF [Related reagents](#)

RRID AB_2271671

Specificity **Sheep anti Human transferrin antibody** recognizes human transferrin, an approximately 80 kDa blood plasma glycoprotein synthesised by the liver, which contains two specific high affinity iron (Fe³⁺) binding sites and is responsible for the transport and supply of an exchangeable pool of iron, through binding to cell surface transferrin receptors.

Transferrin is the primary blood iron transport protein and under normal conditions, approximately one-third of total blood transferrin contains bound iron. Measurement of blood transferrin levels can be used as an indicator for blood iron-carrying capacity and abnormalities of iron metabolism such as anaemia, iron overload and haemochromatosis.

Sheep anti Human transferrin antibody shows minimal cross reactivity with related serum proteins.

References

1. Rouault, T.A. (2003) How mammals acquire and distribute iron needed for oxygen-based metabolism. [PLoS Biol. 1 \(3\): E79.](#)
2. Olkhov, R. and Shaw, A.M. (2014) Growth kinetics of gold nanoparticles on silica/graphene surfaces for multiplex biological immunoassays. [RSC Adv., 2014,4, 31678-84.](#)

Further Reading

1. Giannetti, A.M. *et al.* (2003) Mechanism for multiple ligand recognition by the human transferrin receptor. [PLoS Biol. 1 \(3\): E51.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: 10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Donkey Anti Sheep IgG (STAR88...) [DyLight@488](#), [HRP](#)Rabbit Anti Sheep IgG (H/L) (5184-2304...) [Biotin](#)

North & South Tel: +1 800 265 7376

Worldwide Tel: +44 (0)1865 852 700

Europe Tel: +49 (0) 89 8090 95 21

America Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

To
find a

batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M382558:210513'

Printed on 22 Sep 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)