

Datasheet: 913005

**BATCH NUMBER 155063**

<b>Description:</b>	MOUSE ANTI HUMAN IgA1 ALPHA 1 CHAIN:HRP
<b>Specificity:</b>	IgA1
<b>Format:</b>	HRP
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	B3506B4
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
<b>Preparation</b>	Purified IgG prepared from ascites
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	50% glycerol.
<b>Immunogen</b>	Human IgA1 heavy chains.
<b>External Database Links</b>	<b>UniProt:</b>

[P01876](#)   [Related reagents](#)

**Entrez Gene:**

[3493](#)   IGHA1   [Related reagents](#)

---

**RRID**   AB\_609681

---

**Specificity**   **Mouse anti Human IgA1 antibody, clone B3506B4** recognizes the Fc portion of the heavy chain of human IgA1 as demonstrated by ELISA.

---

**References**

1. Zhao, Y.F. *et al.* (2017) Pathogenic role of glycan-specific IgG antibodies in IgA nephropathy. [BMC Nephrol. 18 \(1\): 301.](#)
2. Li, G. *et al.* (2019) Serum Levels of Joining Chain-Containing IgA1 Are Not Elevated in Patients with IgA Nephropathy. [Dis Markers. 2019: 9802839.](#)

---

**Storage**

Store at -20°C only.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**   Guaranteed until date of expiry. Please see product label.

---

**Health And Safety Information**   Material Safety Datasheet documentation #10106 available at: <https://www.bio-rad-antibodies.com/SDS/913005>  
10106

---

**Regulatory**   For research purposes only

---

## Related Products

### Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

**North & South**   Tel: +1 800 265 7376

**America**   Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M350243:190307'

Printed on 19 Jan 2024