

Datasheet: 913005

BATCH NUMBER 157329

Description:	MOUSE ANTI HUMAN IgA1 ALPHA 1 CHAIN:HRP
Specificity:	IgA1
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	B3506B4
Isotype:	IgG1
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Preparation	Purified IgG prepared from ascites
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	50% glycerol.
Immunogen	Human IgA1 heavy chains.
External Database Links	UniProt:

[P01876](#) [Related reagents](#)

Entrez Gene:

[3493](#) IGHA1 [Related reagents](#)

RRID AB_609681

Specificity **Mouse anti Human IgA1 antibody, clone B3506B4** recognizes the Fc portion of the heavy chain of human IgA1 as demonstrated by ELISA.

References 1. Zhao, Y.F. *et al.* (2017) Pathogenic role of glycan-specific IgG antibodies in IgA nephropathy. [BMC Nephrol. 18 \(1\): 301.](#)
2. Li, G. *et al.* (2019) Serum Levels of Joining Chain-Containing IgA1 Are Not Elevated in Patients with IgA Nephropathy. [Dis Markers. 2019: 9802839.](#)

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10106 available at:
<https://www.bio-rad-antibodies.com/SDS/913005>
10106

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

 Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

 Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

 Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M350243:190307'

Printed on 19 Jan 2024