

Datasheet: 913005

BATCH NUMBER 167035

Description:	MOUSE ANTI HUMAN IgA1 ALPHA 1 CHAIN:HRP
Specificity:	IgA1
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	B3506B4
Isotype:	IgG1
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Preparation	Purified IgG prepared from ascites
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	50% glycerol.
Immunogen	Human IgA1 heavy chains.
External Database Links	UniProt:

[P01876](#) [Related reagents](#)

Entrez Gene:

[3493](#) IGHA1 [Related reagents](#)

RRID AB_609681

Specificity **Mouse anti Human IgA1 antibody, clone B3506B4** recognizes the Fc portion of the heavy chain of human IgA1 as demonstrated by ELISA.

References

1. Zhao, Y.F. *et al.* (2017) Pathogenic role of glycan-specific IgG antibodies in IgA nephropathy. [BMC Nephrol. 18 \(1\): 301.](#)
2. Li, G. *et al.* (2019) Serum Levels of Joining Chain-Containing IgA1 Are Not Elevated in Patients with IgA Nephropathy. [Dis Markers. 2019: 9802839.](#)

Storage

Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information

Material Safety Datasheet documentation #10106 available at:
<https://www.bio-rad-antibodies.com/SDS/913005>
10106

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M401968:220718'

Printed on 19 Jan 2024