

Datasheet: 919001

**BATCH NUMBER 168890**

<b>Description:</b>	MOUSE ANTI HUMAN IgG4
<b>Specificity:</b>	IgG4
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	HP6023
<b>Isotype:</b>	IgG3
<b>Quantity:</b>	0.5 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA	▪			
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared from tissue culture supernatant
<b>Buffer Solution</b>	Borate buffered saline
<b>Preservative Stabilisers</b>	None present
<b>Approx. Protein Concentrations</b>	IgG concentration 0.5 mg/ml
<b>Immunogen</b>	Pool of 5-6 individually purified Human IgG4 myeloma proteins

<b>External Database Links</b>	<b>UniProt:</b> <a href="#">P01861</a> <a href="#">Related reagents</a>
	<b>Entrez Gene:</b> <a href="#">3503</a> IGHG4 <a href="#">Related reagents</a>
<b>RRID</b>	AB_2565207
<b>Specificity</b>	<p><b>Mouse anti Human IgG4 antibody, (clone HP6023)</b> is a monoclonal antibody that recognizes human immunoglobulin G subclass 4 (IgG4). This antibody has been reported to work in the following applications; ELISA, Immunohistology - Frozen.</p> <p>IgG4 is the least abundant human IgG subclass in serum and has been studied in the context of chronic infection and immune mediated disease. Mouse anti Human IgG4 has been used to detect IgG4 antibodies in immunological studies, including investigations of parasitic infection and autoimmune disease.</p>
<b>References</b>	<ol style="list-style-type: none"> <li>Mutapi, F. <i>et al.</i> (2002) The effect of treatment on the age-antibody relationship in children infected with <i>Schistosoma mansoni</i> and <i>Schistosoma haematobium</i>. <a href="#">Mem Inst Oswaldo Cruz. 97 Suppl 1: 173-80.</a></li> <li>Martinsson, K. <i>et al.</i> (2017) Immunoglobulin (Ig)G1 and IgG4 anti-cyclic citrullinated peptide (CCP) associate with shared epitope, whereas IgG2 anti-CCP associates with smoking in patients with recent-onset rheumatoid arthritis (the Swedish TIRA project). <a href="#">Clin Exp Immunol. 188 (1): 53-62.</a></li> </ol>
<b>Storage</b>	<p>Store at -20°C only.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10123 available at: <a href="https://www.bio-rad-antibodies.com/SDS/919001">https://www.bio-rad-antibodies.com/SDS/919001</a>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG IgA IgM (STAR87...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (STAR70...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (STAR77...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (STAR76...)	<a href="#">RPE</a>
Rabbit Anti Mouse IgG (STAR12...)	<a href="#">RPE</a>

Rabbit Anti Mouse IgG (STAR13...) [HRP](#)  
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)  
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)  
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),  
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),  
[FITC](#), [HRP](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M401830:220718'

**Printed on 07 May 2026**

---

© 2026 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)