

Datasheet: 9220-1995

BATCH NUMBER 157274

Description:	NATIVE PIG TRYPSIN
Name:	TRYPSIN
Format:	Semi Pure
Product Type:	Purified Protein
Quantity:	500 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Pig
Product Form	Partially purified enzyme from porcine pancreas - liquid
Preservative Stabilisers	None present
Approx. Protein Concentrations	Total protein concentration 100 mg/ml
External Database Links	UniProt: P00761 Related reagents
Product Information	Partially purified porcine pancreatic trypsin, a serine protease that shows preferential cleavage at Arg- -Xaa, and Lys- -Xaa sites.
Activity	1500 U/mg solid

1 unit causes an increase in absorbance of 0.0003 per minute at 25°C and 253 nm, pH 7.6 on the substrate benzoyl-L arginine ester (BAEE).

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10078 available at:
10078: <https://www.bio-rad-antibodies.com/uploads/MSDS/10078.pdf>

Regulatory For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M363572:200528'

Printed on 22 Mar 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)