

Datasheet: 9400-0039

Description:	RABBIT ANTI BOVINE UBIQUITIN
Specificity:	UBIQUITIN
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA			▪	
Western Blotting	▪			
Immunodiffusion	▪			
Immunofluorescence			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Bovine
Product Form	Ig fraction - lyophilised
Reconstitution	Reconstitute with 1 ml of sterile distilled water. For long term storage the addition of a preservative is recommended.
Preparation	Ig fraction prepared by ammonium sulphate precipitation
Buffer Solution	Ammonium bicarbonate
Preservative Stabilisers	None present

Immunogen Purified bovine ubiquitin.

External Database

Links

UniProt:

[P62990](#)

[Related reagents](#)

RRID

AB_906164

Specificity

Rabbit anti Bovine ubiquitin antibody recognizes bovine ubiquitin.

Storage

Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of reconstitution.

Health And Safety Information

Material Safety Datasheet documentation #10122 available at:
10122: <https://www.bio-rad-antibodies.com/uploads/MSDS/10122.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)

[FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...)

[HRP](#)

Goat Anti Rabbit IgG (Fc) (STAR121...)

[Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...)

[RPE](#)

Sheep Anti Rabbit IgG (STAR36...)

[DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

From March 15, 2021, we will no longer supply printed datasheets with our products.
Look out for updates on how to access your digital version at bio-rad-antibodies.com

'M350792:190314'

Printed on 10 Feb 2021