

## Datasheet: AAI19AB

<b>Description:</b>	SHEEP ANTI BOVINE IgM:Alk. Phos.
<b>Specificity:</b>	IgM
<b>Format:</b>	Alk. Phos.
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.5 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/1000 - 1/10000
Western Blotting	▪			1/5000 - 1/50000

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using the appropriate negative/positive controls.

<b>Target Species</b>	Bovine
<b>Product Form</b>	Purified IgG conjugated to Alkaline Phosphatase - liquid
<b>Antiserum Preparation</b>	Antisera to bovine IgM were raised by repeated immunisation of sheep with highly purified antigen. Purified IgG prepared by affinity chromatography.
<b>Buffer Solution</b>	50mM HEPES, 0.1M NaCl, 1mM MgCl <sub>2</sub> , 0.1mM ZnCl <sub>2</sub>
<b>Preservative Stabilisers</b>	0.09% Sodium Azide 0.2% Bovine Serum Albumin
<b>Approx. Protein Concentrations</b>	IgG concentration 0.5mg/ml
<b>Immunogen</b>	Purified bovine IgM.

<b>RRID</b>	AB_10718859
<b>Specificity</b>	<b>Sheep anti bovine IgM</b> recognizes bovine IgM and shows no cross-reactivity with other bovine immunoglobulin classes in immunoelectrophoresis. Sheep anti bovine IgM may cross-react with IgM from other species.
<b>References</b>	<ol style="list-style-type: none"> <li>1. Nebl, T. <i>et al.</i> (2002) Proteomic analysis of a detergent-resistant membrane skeleton from neutrophil plasma membranes. <a href="#">J Biol Chem. 277 (45): 43399-409.</a></li> <li>2. Assad, A. <i>et al.</i> (2012) Immunophenotyping and characterization of BNP colostrum revealed pathogenic alloantibodies of IgG1 subclass with specificity to platelets, granulocytes and monocytes of all maturation stages. <a href="#">Vet Immunol Immunopathol. 147: 25-34.</a></li> <li>3. Hamsten, C. <i>et al.</i> (2009) Recombinant surface proteomics as a tool to analyze humoral immune responses in bovines infected by <i>Mycoplasma mycoides</i> subsp. <i>mycoides</i> small colony type. <a href="#">Mol Cell Proteomics. 8: 2544-54.</a></li> <li>4. Mansilla FC <i>et al.</i> (2015) Safety and immunogenicity of a soluble native <i>Neospora caninum</i> tachyzoite-extract vaccine formulated with a soy lecithin/β-glucan adjuvant in pregnant cattle. <a href="#">Vet Immunol Immunopathol. 165 (1-2): 75-80.</a></li> <li>5. Hossain, M.M. <i>et al.</i> (2016) Multiplex Detection of IgG and IgM to Rift Valley Fever Virus Nucleoprotein, Nonstructural Proteins, and Glycoprotein in Ovine and Bovine. <a href="#">Vector Borne Zoonotic Dis. Jul 5. [Epub ahead of print]</a></li> <li>6. Van Meulder, F. <i>et al.</i> (2015) Analysis of the protective immune response following intramuscular vaccination of calves against the intestinal parasite <i>Cooperia oncophora</i>. <a href="#">Int J Parasitol. 45 (9-10): 637-46.</a></li> <li>7. Jankowska, A. <i>et al.</i> (2016) Humoral and cellular immune response to <i>Histophilus somni</i> recombinant heat shock protein 60 kDa in farm animals <a href="#">Veterinárni Medicína. 60 (No. 11): 603-12.</a></li> <li>8. Trotta, M. <i>et al.</i> (2015) Simultaneous immunization of cattle with foot-and-mouth disease (FMD) and live anthrax vaccines do not interfere with FMD booster responses <a href="#">Trials in Vaccinology. 4: 38-42.</a></li> </ol>
<b>Storage</b>	<p>Store at +4°C. DO NOT FREEZE..</p> <p>This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10089 available at: 10089: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10089.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10089.pdf</a>
<b>Regulatory</b>	For research purposes only

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

From March 15, 2021, we will no longer supply printed datasheets with our products.  
Look out for updates on how to access your digital version at [bio-rad-antibodies.com](https://www.bio-rad-antibodies.com)

'M363595:200528'

**Printed on 11 Feb 2021**

---

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)