

Datasheet: AAI24

**BATCH NUMBER 158971**

<b>Description:</b>	GOAT ANTI CAT IgM
<b>Specificity:</b>	IgM
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using the appropriate negative/positive controls.

### Target Species

Cat

### Product Form

Purified IgG - liquid

### Antiserum Preparation

Antisera to feline IgM were raised by repeated immunisation of goat with highly purified antigen. Purified IgG prepared by affinity chromatography using antigen coupled to agarose beads.

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>)

### Approx. Protein

IgG concentration 1.0 mg/ml

## Concentrations

**Immunogen** Purified feline IgM.

**RRID** AB\_323017

**Specificity** **Goat anti Cat IgM antibody** recognizes feline IgM and shows no cross - reactivity with other feline immunoglobulin classes in immunoelectrophoresis.

Goat anti Cat IgM antibody may cross react with IgM from other species.

## References

1. Jiménez, A. *et al.* (2008) Membranous glomerulonephritis in the Iberian lynx (*Lynx pardinus*). [Vet Immunol Immunopathol. 121: 34-43.](#)
2. Norris, C.R. *et al.* (2003) Production of polyclonal antisera against feline immunoglobulin E and its use in an ELISA in cats with experimentally induced asthma. [Vet Immunol Immunopathol. 96: 149-57.](#)
3. Limlenglert, P. *et al.* (2011) Application of Flow Cytometry for Early Diagnosis and Monitoring in Cats with Immune-mediated Hemolytic Anemia [Thai J Vet Med. 41: 185-91.](#)

## Storage

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

## Guarantee

12 months from date of despatch

## Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/AAI2410040>

## Regulatory

For research purposes only

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

**North & South America** Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M363613:200528'

Printed on 12 Sep 2024