

Datasheet: AAI25P

BATCH NUMBER 168087

Description:	GOAT ANTI CAT IgA:HRP
Specificity:	IgA
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/10,000 - 1/100,000
Western Blotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using the appropriate negative/positive controls.

Target Species	Cat
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Antiserum Preparation	Antisera to feline IgA were raised by repeated immunisation of goat with highly purified antigen. Purified IgG prepared by affinity chromatography.
Buffer Solution	Phosphate buffered saline.
Preservative Stabilisers	0.05% ProClin 300 0.2% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml.
Immunogen	Purified feline IgA.

RRID	AB_323014
Specificity	Goat anti Cat IgA antibody recognizes feline IgA and shows no cross - reactivity with other feline immunoglobulin classes in immunoelectrophoresis.
References	<ol style="list-style-type: none"> 1. Stevens, R. <i>et al.</i> (2004) Oral immunization with recombinant <i>Listeria monocytogenes</i> controls virus load after vaginal challenge with feline immunodeficiency virus. J Virol. 78 (15): 8210-8. 2. Jiménez, A. <i>et al.</i> (2008) Membranous glomerulonephritis in the Iberian lynx (<i>Lynx pardinus</i>). Vet Immunol Immunopathol. 121 (1-2): 34-43. 3. Balan, P. <i>et al.</i> (2020) Orally administered ovine serum immunoglobulins modulate dental plaque in cats. Res Vet Sci. 133: 262-8. 4. Gore, A.M. <i>et al.</i> (2021) Supplementation of Diets With Bovine Colostrum Influences Immune and Gut Function in Kittens. Front Vet Sci. 8: 675712. 5. Mohamed Cherif, A. <i>et al.</i> (2022) A cross- sectional study of Leishmania infantum infection in stray cats in Algiers' suburbs, Algeria, and evaluation of serological and molecular tests for its diagnosis. Veterinaria 71 (1): 73-83.
Storage	<p>Store at +4°C. DO NOT FREEZE.</p> <p>Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch.
Acknowledgements	ProClin is a trademark of The Dow Chemical Company ("Dow") or an affiliated company of Dow.
Health And Safety Information	<p>Material Safety Datasheet documentation #20391 available at: https://www.bio-rad-antibodies.com/SDS/AAI25P</p> <p>20391</p>
Regulatory	For research purposes only.

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M428049:240301'

Printed on 12 Sep 2024

