

## Datasheet: AAI48

**BATCH NUMBER 156668**

<b>Description:</b>	GOAT ANTI PIG IgM
<b>Specificity:</b>	IgM
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			1/200 - 1/2000
Immunohistology - Paraffin			▪	
ELISA (1)	▪			1/1000 - 1/30,000
Immunoprecipitation			▪	
Western Blotting	▪			1/1000 - 1/30,000
Immunocytochemistry			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1)For coating plates a 1/100 - 1/500 dilution is recommended.

<b>Target Species</b>	Pig
<b>Product Form</b>	Purified IgG - liquid
<b>Antiserum Preparation</b>	Antisera to porcine IgM were raised by repeated immunisation of goat with highly purified antigen. Purified IgG prepared by affinity chromatography using antigen coupled to agarose beads.
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )

<b>Approx. Protein Concentrations</b>	IgG concentration 1.0 mg/ml
<b>Immunogen</b>	Purified Porcine IgM
<b>Specificity</b>	<b>Goat anti Pig IgM antibody</b> recognizes porcine IgM and shows no cross-reactivity with other porcine immunoglobulin classes in immunoelectrophoresis. This antibody may cross-react with IgM from other species.
<b>References</b>	<ol style="list-style-type: none"> <li>Williams, A.R. <i>et al.</i> (2017) Dietary cinnamaldehyde enhances acquisition of specific antibodies following helminth infection in pigs. <a href="#">Vet Immunol Immunopathol. 189: 43-52.</a></li> <li>Tiurbe, G. <i>et al.</i> (2009) Inhibitory effects of rat bone marrow-derived dendritic cells on naïve and alloantigen-specific CD4+ T cells: a comparison between dendritic cells generated with GM-CSF plus IL-4 and dendritic cells generated with GM-CSF plus IL-10. <a href="#">BMC Res Notes. 2: 12.</a></li> <li>Corsaut, L. <i>et al.</i> (2021) Immunogenicity study of a <i>Streptococcus suis</i>. autogenous vaccine in preparturient sows and evaluation of passive maternal immunity in piglets. <a href="#">BMC Vet Res. 17 (1): 72.</a></li> <li>Forner, R. <i>et al.</i> (2021) Distribution difference of colostrum-derived B and T cells subsets in gilts and sows. <a href="#">PLoS One. 16 (5): e0249366.</a></li> </ol>
<b>Storage</b>	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted.
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/SDS/AAI4810040">https://www.bio-rad-antibodies.com/SDS/AAI4810040</a>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M363672:200528'

Printed on 01 Mar 2024