

Datasheet: AAI50P

BATCH NUMBER 167844

Description:	GOAT ANTI DOG IgG (H+L):HRP
Specificity:	IgG
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			1/4000 - 1/8000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Dog
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid

Antiserum Preparation Purified IgG from goat antisera prepared by affinity chromatography on canine IgG covalently linked to agarose.

Buffer Solution Phosphate buffered saline.

Preservative Stabilisers 50% Glycerol.

Immunogen Canine IgG

Specificity **Goat anti Dog IgG (H+L)** recognizes the heavy and light chains of canine IgG and may react with the light chains of other canine immunoglobulins.

References 1. Rupp, A. *et al.* (2013) Anti-GM2 ganglioside antibodies are a biomarker for acute canine polyradiculoneuritis. [J Peripher Nerv Syst. 18 \(1\): 75-88.](https://doi.org/10.1007/s12035-013-0388-8)

2. Gianesini, G. *et al.* (2023) Immune-Mediated Hemolytic Anemia and Clinically Suspected Acute Pancreatitis in Dogs, a Pilot Study. [Top Companion Anim Med. 56-57: 100821.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

Guarantee Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information Material Safety Datasheet documentation #10106 available at: <https://www.bio-rad-antibodies.com/SDS/AAI50P>
10106

Regulatory For research purposes only.

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M420558:230706'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)