

Datasheet: AAM13B

Description:	RABBIT ANTI MOUSE INTERLEUKIN-1 BETA:Biotin
Specificity:	IL-1 BETA
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.25 - 1.0ug/ml
Western Blotting	▪			0.1 - 0.2ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using the appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgG conjugated to Biotin - lyophilized
Reconstitution	Reconstitute with 0.5 ml sterile PBS containing 0.1% bovine serum albumin Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended.
Antiserum Preparation	Antisera to mouse IL-1 beta were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present

Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml after reconstitution
Immunogen	Recombinant murine IL-1 beta.
External Database Links	<p>UniProt: P10749 Related reagents</p> <p>Entrez Gene: 16176 Il1b Related reagents</p>
RRID	AB_2124617
Specificity	Rabbit anti Mouse Interleukin-1 beta antibody recognizes mouse Interleukin-1 beta (IL-1 beta), a cytokine involved in the inflammatory response.
ELISA	AAM13B may be used in a direct ELISA or as the detection reagent in a sandwich ELISA with AAM13G as the capture reagent and PMP71 as the standard.
Western Blotting	AAM13B may be used under either reducing or non-reducing conditions with PMP71 as the positive control.
References	<ol style="list-style-type: none"> 1. Clausen, B.H. <i>et al.</i> (2008) Interleukin-1beta and tumor necrosis factor-alpha are expressed by different subsets of microglia and macrophages after ischemic stroke in mice. J Neuroinflammation. 5: 46. 2. Ruan, L. <i>et al.</i> (2009) Amyloid deposition and inflammation in APP^{swe}/PS1^{dE9} mouse model of Alzheimer's disease. Curr Alzheimer Res. 6: 531-40.
Storage	<p>Prior to reconstitution store at -20°C. After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	6 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10162 available at: 10162: https://www.bio-rad-antibodies.com/uploads/MSDS/10162.pdf
Regulatory	For research purposes only

Related Products

Recommended Positive Controls

[RECOMBINANT MOUSE INTERLEUKIN-1 BETA \(PMP71\)](#)

Recommended Useful Reagents

RABBIT ANTI MOUSE INTERLEUKIN-1 BETA (AAM13G)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M391248:211008'

Printed on 08 Oct 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)