

Datasheet: AAM17G

Description:	RABBIT ANTI MOUSE FGF 9
Specificity:	FGF 9
Other names:	FIBROBLAST GROWTH FACTOR 9
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.5 ug/ml
Western Blotting	▪			0.1ug/ml - 0.2ug/ml
Functional Assays	▪			25ng/ml - 60ng/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using the appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgG - lyophilised
Reconstitution	Reconstitute with 0.1 ml distilled water
Antiserum Preparation	Antisera to murine FGF 9 were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared by antigen affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present.
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml after reconstitution

Immunogen	Recombinant murine FGF-9.
External Database Links	<p>UniProt: P54130 Related reagents</p> <p>Entrez Gene: 14180 Fgf9 Related reagents</p>
Synonyms	Fgf-9
RRID	AB_2278293
Specificity	Rabbit anti Mouse FGF-9 antibody recognizes mouse Fibroblast Growth Factor 9 (FGF 9). Rabbit anti Mouse FGF-9 antibody neutralizes murine FGF-9 activity with an ND ₅₀ for mL-6 (1.5 ngml ⁻¹ of between 0.025 and 0.06µg ml ⁻¹ .
Storage	<p>Prior to reconstitution store at -20°C. After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10162 available at: 10162: https://www.bio-rad-antibodies.com/uploads/MSDS/10162.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
---	------------------	---	---------------	---

'M373007:200812'

Printed on 12 Aug 2020