

## Datasheet: AAM35

<b>Description:</b>	RABBIT ANTI MOUSE INTERLEUKIN-2
<b>Specificity:</b>	IL-2
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

## Product Details

**RRID** AB\_324469

**Applications** This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.5 - 2.0ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1 - 0.2ug/ml
Functional Assays	▪			0.04 - 0.05ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using the appropriate negative/positive controls.

**Target Species** Mouse

**Product Form** Purified IgG - lyophilised

**Reconstitution** Reconstitute with 0.1 ml distilled water  
Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended.  
N.B. For functional studies do not add sodium azide

**Antiserum Preparation** Antisera to mouse IL-2 were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

**Buffer Solution** Phosphate buffered saline

**Preservative Stabilisers** None present

<b>Approx. Protein Concentrations</b>	IgG concentration 1.0 mg/ml after reconstitution
<b>Immunogen</b>	<a href="#">Recombinant mouse IL-2</a> (PMP38)
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P04351</a>   <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">16183</a>   IL2   <a href="#">Related reagents</a></p>
<b>Synonyms</b>	IL-2
<b>Specificity</b>	<p><b>Rabbit anti Mouse Interleukin-2 antibody</b> recognizes mouse Interleukin-2 (IL-2), a potent lymphoid cell growth factor which primarily acts on T lymphocytes, in response to antigenic or mitogenic stimulation.</p> <p>IL-2 stimulates the growth/differentiation of monocytes, macrophages, B lymphocytes, NK cells, and lymphokine activated killer cells which are used in experimental cancer therapy.</p>
<b>ELISA</b>	AAM35 may be used in an indirect ELISA or as a capture antibody in a sandwich ELISA with <a href="#">AAM35B</a> as the detection reagent and <a href="#">PMP38</a> as the standard.
<b>References</b>	1. Smith, K.A. (1998) Interleukin-2: inception, impact, and implications. <a href="#">Science, 240: 1169-1176.</a>
<b>Storage</b>	<p>Prior to reconstitution store at +4°C. Following reconstitution store at -20°C.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	12 months from date of reconstitution.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10162 available at: 10162: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10162.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10162.pdf</a>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...)   [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...)   [RPE](#)
- Sheep Anti Rabbit IgG (2AB02...)   [Biotin](#)
- Sheep Anti Rabbit IgG (STAR36...)   [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

### Recommended Useful Reagents

TidyBlot WESTERN BLOT DETECTION REAGENT:HRP (STAR209P)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

'M341200:190109'

**Printed on 11 Oct 2019**

---

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)