

Datasheet: AAM51

Description:	RABBIT ANTI MOUSE VEGF
Specificity:	VEGF
Other names:	VPF
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			2.5ug/ml
ELISA	▪			0.5ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1ug/ml - 0.2ug/ml
Immunofluorescence	▪			2.5ug/ml
Functional Assays	▪			1.0ug/ml - 2.0ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. 1mM EDTA buffer pH 8.0 is recommended for this purpose.

Target Species	Mouse
Product Form	Purified IgG - lyophilised
Reconstitution	<p>Reconstitute with 0.1ml distilled water</p> <p>Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended.</p> <p>N.B. For functional studies do not add sodium azide</p>
Antiserum Preparation	Antisera to mouse VEGF were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative	None present

Stabilisers

Carrier Free Yes

Approx. Protein Concentrations IgG concentration 1.0 mg/ml after reconstitution.

Immunogen [Recombinant mouse VEGF](#)

External Database Links

UniProt:

[Q00731](#) [Related reagents](#)

Entrez Gene:

[22339](#) Vegfa [Related reagents](#)

Synonyms Vegf

RRID AB_2288247

Specificity **Rabbit anti Mouse VEGF antibody** recognizes mouse Vascular Endothelial Growth Factor (VEGF), a protein active in angiogenesis, vasculogenesis and endothelial cell growth.

Rabbit anti Mouse VEGF antibody has not been cross adsorbed and may react with homologous proteins from related species.

ELISA This antibody may be used in indirect ELISA or as a capture reagent in a sandwich ELISA with [AAM51B](#) as a detection reagent and [PMP67](#) as the standard.

Histology Positive Control Tissue Kidney

References

1. Nahrendorf, M. *et al.* (2007) The healing myocardium sequentially mobilizes two monocyte subsets with divergent and complementary functions. [J Exp Med. 204: 3037-47.](#)
 2. Araújo, W.F. *et al.* (2015) Renin-angiotensin system (RAS) blockade attenuates growth and metastatic potential of renal cell carcinoma in mice. [Urol Oncol. 33 \(9\): 389.e1-7.](#)
 3. Lin, E.Y. *et al.* (2007) Vascular endothelial growth factor restores delayed tumor progression in tumors depleted of macrophages. [Mol Oncol. 1: 288-302.](#)
 4. HogenEsch H *et al.* (2016) Angiogenesis in the skin of SHARPIN-deficient mice with chronic proliferative dermatitis. [Exp Mol Pathol. Oct 26. pii: S0014-4800\(16\)30084-3. \[Epub ahead of print\]](#)
 5. Hogenesch, H. *et al.* (2016) Angiogenesis in the skin of SHARPIN-deficient mice with chronic proliferative dermatitis. [Exp Mol Pathol. 101 \(3\): 303-7.](#)
-

Storage Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.

This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 6 months from date of despatch

Health And Safety Material Safety Datasheet documentation #10162 available at:

Information 10162: <https://www.bio-rad-antibodies.com/uploads/MSDS/10162.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Sheep Anti Rabbit IgG (STAR36...) [DyLight@488](#), [DyLight@680](#), [DyLight@800](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)
[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025C\)](#)
[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

'M372913:200818'

Printed on 18 Aug 2020

© 2020 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)