

## Datasheet: AAR16G

**BATCH NUMBER 156560**

<b>Description:</b>	RABBIT ANTI RAT INTERLEUKIN-4
<b>Specificity:</b>	IL-4
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.5 ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1ug/ml - 0.2ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using the appropriate negative/positive controls.

#### Target Species

Rat

#### Product Form

Purified IgG - lyophilised

#### Reconstitution

Reconstitute with 0.1ml distilled water.

Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended

#### Antiserum Preparation

Antisera to rat IL-4 were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.

#### Buffer Solution

Phosphate buffered saline

<b>Preservative Stabilisers</b>	None present.
<b>Carrier Free</b>	Yes
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0 mg/ml after reconstitution.
<b>Immunogen</b>	Recombinant rat IL-4.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P20096</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">287287</a>    IL4    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	IL-4
<b>RRID</b>	AB_323569
<b>Specificity</b>	<b>Rabbit anti Rat Interleukin-4 antibody</b> specifically recognizes rat interleukin-4 (IL-4), also known as B-cell stimulatory factor 1, a pleiotropic cytokine which interacts with multiple cell lineages, including T cells, fibroblasts and haematopoietic cells, and is one of the most important factors involved in B cell IgG1 and IgE class switching. IL-4 acts as a regulator of low affinity Fc receptor for IgE (CD23) expression on both lymphocytes and monocytes.
<b>Storage</b>	<p>Prior to reconstitution store at -20°C.  After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10162 available at: <a href="https://www.bio-rad-antibodies.com/SDS/AAR16G">https://www.bio-rad-antibodies.com/SDS/AAR16G</a> 10162
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...)    [FITC](#)  
Goat Anti Rabbit IgG (Fc) (STAR121...)    [Biotin](#), [FITC](#), [HRP](#)  
Sheep Anti Rabbit IgG (STAR35...)    [RPE](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

## Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M372435:200706'

**Printed on 18 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)