

Datasheet: AAR17B

BATCH NUMBER 160823

Description:	RABBIT ANTI RAT INTERLEUKIN-6:Biotin
Specificity:	IL-6
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1.0 ug/ml
Western Blotting	▪			0.1ug/ml - 0.2ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rat
Product Form	Purified IgG conjugated to Biotin - lyophilised
Reconstitution	Reconstitute with 0.5 ml sterile PBS containing 0.1% bovine serum albumin Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended.
Antiserum Preparation	Antisera to rat IL-6 were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative	None Present

Stabilisers

Carrier Free Yes

Approx. Protein Concentrations IgG concentration 0.1 mg/ml after reconstitution

Immunogen [Recombinant Rat IL-6](#) (PRP21)

External Database Links

UniProt:
[P20607](#) [Related reagents](#)

Entrez Gene:
[24498](#) IL6 [Related reagents](#)

Synonyms IL-6

RRID AB_2249219

Specificity **Rabbit anti Rat Interleukin-6 antibody** recognizes rat interleukin-6.

IL-6 is a pleiotropic cytokine involved in the regulation of immune and inflammatory responses.

ELISA Anti Rat IL-6: biotin conjugate may be used in a direct ELISA with [recombinant rat IL-6](#) (PRP21) as the standard.

Storage Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.

This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10162 available at:
<https://www.bio-rad-antibodies.com/SDS/AAR17B>
10162

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M372436:200706'

