

## Datasheet: AAR42

**BATCH NUMBER 159852**

<b>Description:</b>	RABBIT ANTI RAT VEGF
<b>Specificity:</b>	VEGF
<b>Other names:</b>	VPF
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			0.5 - 2.0ug/ml
Immunoprecipitation			■	
Western Blotting	■			0.1 - 0.2ug/ml
Functional Assays	■			0.05 - 0.1ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Rat
<b>Product Form</b>	Purified IgG - lyophilised
<b>Reconstitution</b>	<p>Reconstitute with 0.1ml distilled water</p> <p>Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended.</p>
<b>Antiserum Preparation</b>	Antiserum to rat VEGF was raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared by affinity chromatography.

<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	None present
<b>Carrier Free</b>	Yes
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0mg/ml after reconstitution
<b>Immunogen</b>	Recombinant rat VEGF.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P16612</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">83785</a>    Vegfa    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	Vegf
<b>RRID</b>	AB_2212685
<b>Specificity</b>	<p><b>Rabbit anti Rat VEGF antibody</b> specifically recognizes rat VEGF (vascular endothelial growth factor), a heparin-binding cytokine and member of the PDGF/VEGF growth factor family, which acts as a vasodilator, regulator of vasculogenesis, and is a potent stimulator of angiogenesis of both normal and cancerous cells.</p> <p>All VEGF family members act through binding to the tyrosine kinase receptors VEGFR-1 (Flt-1) and VEGFR-2 (KDR/CD309), and the functional forms of VEGF, which exist as disulphide-linked homodimers, are classified according to the number of amino acid residues within a constitutive polypeptide chain e.g VEGF<sub>165</sub>.</p>
<b>ELISA</b>	Purified rabbit anti rat VEGF may be used in an indirect ELISA or as the capture reagent in a sandwich ELISA with our <a href="#">biotinylated VEGF antibody</a> as the detection reagent.
<b>Storage</b>	<p>Prior to reconstitution store at -20°C.  After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	6 months from date of despatch
<b>Health And Safety Information</b>	<p>Material Safety Datasheet documentation #10162 available at:  <a href="https://www.bio-rad-antibodies.com/SDS/AAR42">https://www.bio-rad-antibodies.com/SDS/AAR42</a>  10162</p>

## Related Products

### Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

### Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M390870:211006'

Printed on 18 Jan 2024