

Datasheet: AAR42B

Description:	RABBIT ANTI RAT VEGF:Biotin
Specificity:	VEGF
Other names:	VPF
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

RRID AB_2212684

Applications This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			0.5 - 2.0ug/ml
Immunoprecipitation			■	
Western Blotting	■			0.1 - 0.2ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Rat

Product Form Purified IgG conjugated to Biotin - lyophilised

Reconstitution Reconstitute with 0.5ml sterile PBS containing 0.1% bovine serum albumin. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended.

Antiserum Preparation Antiserum to rat VEGF was raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared by affinity chromatography.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers None present

Approx. Protein Concentrations	IgG concentration 0.1mg/ml after reconstitution
Immunogen	Recombinant rat VEGF.
External Database Links	<p>UniProt: P16612 Related reagents</p> <p>Entrez Gene: 83785 Vegfa Related reagents</p>
Synonyms	Vegf
Specificity	<p>Rabbit anti Rat VEGF antibody specifically recognizes rat VEGF (vascular endothelial growth factor), a heparin-binding cytokine and member of the PDGF/VEGF growth factor family, which acts as a vasodilator, regulator of vasculogenesis, and is a potent stimulator of angiogenesis of both normal and cancerous cells.</p> <p>All VEGF family members act through binding to the tyrosine kinase receptors VEGFR-1 (Flt-1) and VEGFR-2 (KDR/CD309), and the functional forms of VEGF, which exist as disulphide-linked homodimers, are classified according to the number of amino acid residues within a constitutive polypeptide chain e.g VEGF₁₆₅.</p>
ELISA	Biotinylated rabbit anti rat VEGF may be used in a direct ELISA or as the detection reagent in a sandwich ELISA with our purified rabbit anti rat VEGF antibody as the capture reagent.
Storage	<p>Prior to reconstitution store at +4°C. After reconstitution store at -20°C.</p> <p>Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Shelf Life	12 months from date of reconstitution.
Health And Safety Information	Material Safety Datasheet documentation #10162 available at: 10162: https://www.bio-rad-antibodies.com/uploads/MSDS/10162.pdf
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

'M338970:181217'

Printed on 02 Jan 2019