

Datasheet: AHP1001

BATCH NUMBER 158463

Description:	RABBIT ANTI TSC1
Specificity:	TSC1
Other names:	TUBEROUS SCLEROSIS COMPLEX 1
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			1ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Species Cross Reactivity

Reacts with: Rat, Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Antiserum Preparation Antisera to human TSC1 were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.02% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	A 15 amino acid peptide located near the centre of human TSC1.
External Database Links	<p>UniProt: Q92574 Related reagents</p> <p>Entrez Gene: 7248 TSC1 Related reagents</p>
Synonyms	KIAA0243, TSC
RRID	AB_609656
Specificity	<p>Rabbit anti Human TSC1 antibody recognizes an epitope within the internal region of the 130kDa tumour suppressor protein, Tumour sclerosis complex 1 (TSC1). Mutation within the TSC1 gene is responsible for Tuberous sclerosis (TSC), an autosomal dominant genetic disorder associated with the development of benign tumors (hamartomas), a condition also induced by mutation of the Tumour sclerosis complex 2 (TSC2) gene.</p> <p>Many proteins have been shown to play a role in the regulation of the serine/threonine kinase known as target of rapamycin (TOR), a central component in a complex signaling pathway which controls cell proliferation and cell cycle progression. Negative regulation of TOR activity occurs following the over-expression of the TSC1 and TSC2 gene products, hamartin and tuberin. These act by suppressing the phosphorylation of eukaryotic translation initiation factor 4E binding protein 1 (4EBP1) and ribosomal protein S6 kinase (S6K), downstream targets of TOR involved in mRNA translation. Furthermore the TOR activator protein Rheb has been identified as a target of TSC1-TSC2 complexes.</p> <p>Rabbit anti Human TSC1 antibody has been reported as suitable for use in IHC on mouse EL4 cell cytospins at 2µg/ml.</p>
Western Blotting	AHP1001 detects a band of approximately 130kDa in mouse EL4 cell lysates.
Further Reading	<ol style="list-style-type: none"> Inoki, K. <i>et al.</i> (2005) Signaling by target of rapamycin proteins in cell growth control. Microbiol Mol Biol Rev. 69 (1): 79-100. Tee, A.R. <i>et al.</i> (2002) Tuberous sclerosis complex-1 and -2 gene products function together to inhibit mammalian target of rapamycin (mTOR)-mediated downstream signaling. Proc Natl Acad Sci U S A. 99 (21): 13571-6.
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for

short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP1001 10040
--------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Recommended Useful Reagents

- [ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)
- [TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M386145:210519'

Printed on 18 Jan 2024