

Datasheet: AHP1009

BATCH NUMBER 163910

Description:	RABBIT ANTI BECLIN-1 (N-TERMINAL)
Specificity:	BECLIN-1 (N-TERMINAL)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			0.5ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species Human

Species Cross Reactivity Reacts with: Mouse, Rat
N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form Purified IgG - liquid

Antiserum Preparation Antisera to human Beclin-1 (NT) were raised by repeated immunisations of rabbits with

highly purified antigen. Purified IgG prepared by affinity chromatography.

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative Stabilisers	0.02% Sodium Azide (NaN ₃)
---------------------------------	--

Approx. Protein Concentrations	IgG concentration 1.0mg/ml
---------------------------------------	----------------------------

Immunogen	A 17 amino acid peptide located near human Beclin-1 amino-terminus.
------------------	---

External Database Links	UniProt: Q14457 Related reagents O88597 Related reagents Q91XJ1 Related reagents Entrez Gene: 8678 BECN1 Related reagents 56208 Becn1 Related reagents 114558 Becn1 Related reagents
--------------------------------	---

Synonyms	GT197
-----------------	-------

RRID	AB_609563
-------------	-----------

Specificity	<p>Rabbit anti Human Beclin-1 antibody recognizes an epitope within the N-terminal region of Beclin-1, a ~52 kDa coiled-coil protein which shares structural homology with the yeast autophagy gene <i>apg6/vps30</i>.</p> <p>The ability of cells to undergo bulk degradation of endogenous proteins through an autophagosomic-lysosomal pathway, termed autophagy or Type II cell death, is important in the control of cell growth and development. Beclin-1 is a Bcl-2 interacting mammalian gene which not only mediates autophagy, but also functions as a tumor suppressor, being expressed at reduced levels in human breast carcinoma, and may also play a role in defense against viral infection.</p>
--------------------	---

Western Blotting	AHP1009 detects a band of approximately 52kDa in human 293 cell lysates.
-------------------------	--

Further Reading	<ol style="list-style-type: none">1. Liang, X.H. <i>et al.</i> (1999) Induction of autophagy and inhibition of tumorigenesis by beclin 1. Nature. 402 (6762): 672-6.2. Gozuacik, D. & Kimchi, A. (2004) Autophagy as a cell death and tumor suppressor mechanism. Oncogene. 23 (16): 2891-906.3. Liang, X.H. <i>et al.</i> (1998) Protection against fatal Sindbis virus encephalitis by beclin, a novel Bcl-2-interacting protein. J Virol. 72 (11): 8586-96.
------------------------	--

Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at
----------------	---

-20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP100910040
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)
[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M386223:210519'

Printed on 25 Mar 2023