

Datasheet: AHP1011

BATCH NUMBER 166184

Description:	RABBIT ANTI HUMAN CLUSTERIN
Specificity:	CLUSTERIN
Other names:	APOLIPOPROTEIN J
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			0.5 - 1.0ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species	Human
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to Clusterin were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.
Buffer Solution	Phosphate buffered saline

Preservative Stabilisers	0.02% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	Recombinant human clusterin isoform 1.
External Database Links	<p>UniProt: P10909 Related reagents</p> <p>Entrez Gene: 1191 CLU Related reagents</p>
Synonyms	APOJ, CLI, KUB1
RRID	AB_2276520
Specificity	<p>Rabbit anti Human clusterin antibody recognizes human clusterin, also known as apolipoprotein J, Ku70-binding protein 1 or TRPM-2. Clusterin is a ubiquitous glycoprotein which exists in at least two distinct isoforms; a proapoptotic nuclear form, of ~49 kDa, and a prosurvival secretory form with a molecular mass of 75-80 kDa. Clusterin is implicated in several diverse physiological processes including sperm maturation, cell differentiation, and apoptosis.</p> <p>Clusterin upregulation has been reported in various human malignancies, including bladder, prostate, and breast tumours.</p> <p>Following specific stress induced injury, secretory clusterin acts as a molecular chaperone, scavenging denatured proteins outside cells. Clusterin upregulated in Alzheimer's disease (Miners et al. 2017), retinitis pigmentosa (Jomary et al 1993), glomerulonephritis, atherosclerosis and myocardial infarction (Rosenberg et al. 2002).</p>
Western Blotting	AHP1011 detects a band of approximately 60kDa in human brain tissue lysate.
Further Reading	<ol style="list-style-type: none"> 1. Trougakos, I.P. & Gonos, E.S. (2002) Clusterin/apolipoprotein J in human aging and cancer. Int J Biochem Cell Biol. 34 (11): 1430-48. 2. Calero, M. et al. (2005) Clusterin and Alzheimer's disease. Subcell Biochem. 38: 273-98. 3. Leskov, K.S. et al. (2003) Synthesis and functional analyses of nuclear clusterin, a cell death protein. J Biol Chem. 278 (13): 11590-600. 4. Carreras, I. et al. (2005) Upregulation of clusterin/apolipoprotein J in lactacystin-treated SH-SY5Y cells. J Neurosci Res. 79 (4): 495-502.
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP1011 10040
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)
[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M381669:210512'

Printed on 18 Jan 2024