

Datasheet: AHP1032B

BATCH NUMBER 165067

Description:	RABBIT ANTI HUMAN INTERLEUKIN-2:Biotin
Specificity:	IL-2
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.25 - 1.0ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1 - 0.2ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Product Form

Purified IgG conjugated to Biotin - lyophilized

Reconstitution

Reconstitute with 0.5 ml sterile PBS containing 0.1% Bovine Serum Albumin. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended.

Antiserum Preparation

Antisera to human IL-2 were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers	None present
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml after reconstitution.
Immunogen	Recombinant human IL-2
External Database Links	<p>UniProt: P60568 Related reagents</p> <p>Entrez Gene: 3558 IL2 Related reagents</p>
RRID	AB_2125582
Specificity	<p>Rabbit anti Human Interleukin-2 antibody recognizes human IL-2 (Interleukin-2), otherwise known as T-cell growth factor, a 133 amino acid glycoprotein and potent immunoregulatory cytokine produced by antigen or mitogen stimulated T cells.</p> <p>IL-2 is a multi-functional cytokine which signals through the IL-2 receptor complex (IL-2R), and is principally involved in the stimulation of the T cell proliferative response, promoting T cell clonal expansion and the expansion of T memory cells expressing low levels of high-affinity IL-2R.</p> <p>IL-2 also stimulates the growth/differentiation of B cells, monocytes, NK cells, macrophages and lymphokine activated killer cells.</p>
ELISA	This biotinylated human IL-2 antibody may be used in a direct ELISA or as the detection reagent in a sandwich ELISA with our purified human IL-2 antibody (AHP1032) and recombinant human IL-2 (PHP042) as the standard.
Further Reading	<ol style="list-style-type: none"> 1. Fujita, T. <i>et al.</i> (1983) Structure of the human interleukin 2 gene. Proc Natl Acad Sci U S A. 80 (24): 7437-41. 2. Nelson, B.H. (2004) IL-2, regulatory T cells, and tolerance. J Immunol. 172 (7): 3983-8.
Storage	<p>Prior to reconstitution store at -20°C. After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety	Material Safety Datasheet documentation #10294 available at:

Information <https://www.bio-rad-antibodies.com/SDS/AHP1032B>
10294

Regulatory For research purposes only

Related Products

Recommended Positive Controls

[RECOMBINANT HUMAN INTERLEUKIN-2 \(PHP042A\)](#)

Recommended Useful Reagents

[RABBIT ANTI HUMAN INTERLEUKIN-2 \(AHP1032\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M402198:220718'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)