

Datasheet: AHP1044

Description:	RABBIT ANTI 14-3-3 BETA
Specificity:	14-3-3 BETA
Format:	Serum
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			1/600
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting	▪			1/3000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections.

Target Species

Sheep

Species Cross Reactivity

Reacts with: Mouse, Rabbit, Human, Rat, Bovine

Does not react with: Amphibia, Fungal, Fish, Invertebrates, Reptile, Yeast

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Serum - liquid

Antiserum Preparation Antisera to sheep 14-3-3 beta were raised by repeated immunisations of rabbit with highly

purified antigen.

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)

Immunogen

Synthetic peptide corresponding to acetylated N-terminal sequence of sheep 14-3-3 beta.

Peptide sequence: Ac.TMDKSELVC

External Database Links

UniProt:

[P68251](#)

[Related reagents](#)

RRID

AB_2217621

Specificity

Rabbit anti 14-3-3 β antibody recognizes the acetylated N-terminal of 14-3-3 β. It is a member of the 14-3-3 family which consists of 30 kDa proteins that are involved in multiple protein kinase signalling pathways, regulation of cell cycle progression, cytoskeletal structure, transcription, intracellular trafficking and targeting. Protein interactions with 14-3-3 show distinct preference for its different isoforms and are regulated by phosphorylation of both 14-3-3 and the bound protein.

14-3-3 β plays a role in the pathogenesis of systemic lupus erythematosus and is involved in negative regulation of p90 ribosomal S6 kinase activity and tuberous sclerosis tumour suppressor proteins. It also binds IGF-1 receptor, RAF1 and CDC25 phosphatases.

Rabbit anti 14-3-3 β antibody may not react with recombinant proteins that are not N-acetylated.

Histology Positive Control Tissue

Normal brain tissue

Western Blotting

AHP1044 detects a band of approximately 30kDa in HEK293 cell lysates

References

1. Martin, H. *et al.* (1993) Antibodies against major brain isoforms of 14-3-3 protein. An antibody specific for the N-acetylated amino-terminus of a protein. [FEBS Lett. 331: 296-303.](#)
2. Chen, H.K. *et al.* (2003) Interaction of Akt-phosphorylated ataxin-1 with 14-3-3 mediates neurodegeneration in spinocerebellar ataxia type 1. [Cell. 113: 457-68.](#)
3. Wang, Z. *et al.* (2011) The prognostic value of 14-3-3 isoforms in vulvar squamous cell carcinoma cases: 14-3-3β and ε are independent prognostic factors for these tumors. [PLoS One. 6: e24843.](#)
4. Baxter, H.C. *et al.* (2002) Immunolocalisation of 14-3-3 isoforms in normal and scrapie-infected murine brain. [Neuroscience. 109: 5-14.](#)
5. Chaudhri, M. *et al.* (2003) Mammalian and yeast 14-3-3 isoforms form distinct patterns of dimers in vivo. [Biochem Biophys Res Commun. 300: 679-85.](#)
6. Taurino, F. *et al.* (2014) Short-term Type-1 diabetes differentially modulates 14-3-3 proteins in rat brain and liver. [Eur J Clin Invest. 44: 350-8.](#)

Further Reading 1. Aitken, A. (2006) 14-3-3 proteins: a historic overview. [Semin Cancer Biol. 16 \(3\): 162-72.](#)

Storage Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 18 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10081 available at:
10081: <https://www.bio-rad-antibodies.com/uploads/MSDS/10081.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

From March 15, 2021, we will no longer supply printed datasheets with our products.

Look out for updates on how to access your digital version at bio-rad-antibodies.com

'M352390:190325'

Printed on 09 Feb 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)