

Datasheet: AHP1074

BATCH NUMBER 165820

Description:	RABBIT ANTI DYKDDDDK TAG
Specificity:	DYKDDDDK TAG
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Synthetic Peptide

Product Form Purified IgG - liquid

Antiserum Preparation Antisera to peptide DYKDDDDK were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers <0.1% Sodium Azide (NaN₃)

Approx. Protein Concentrations IgG concentration 0.5mg/ml

Immunogen	Keyhole limpet haemocyanin conjugated Enterokinase cleavage site peptide DYKDDDDK. Residues of cysteine and glycine were used to facilitate coupling at the C-terminal end.
RRID	AB_770097
Specificity	<p>Rabbit anti DYKDDDDK Tag antibody recognizes the hydrophilic epitope tagging octapeptide sequence DYKDDDDK (Asp-Tyr-Lys-Asp-Asp-Asp-Lys), when fused to either the amino- or carboxyl- terminus of target proteins.</p> <p>Rabbit anti DYKDDDDK Tag antibody has been tested against both the immunogen and recombinant proteins containing the DYKDDDDK sequence in ELISA and Western blotting and shows greater binding affinity and sensitivity to fusion proteins containing the DYKDDDDK sequence than the clones M1, M2 and M5, in the majority of assays.</p> <p>In Western blotting of bacterial extracts Rabbit anti DYKDDDDK Tag antibody does not cross-react with any endogenous proteins.</p>
References	1. Keller, B. <i>et al.</i> (2016) Early onset combined immunodeficiency and autoimmunity in patients with loss-of-function mutation in LAT. J Exp Med. May 30. pii: jem.20151110. [Epub ahead of print]
Further Reading	1. Hopp, T.P. <i>et al.</i> (1988) A short polypeptide marker sequence useful for recombinant protein identification and purification. Bio/Technology 6: 1204-10.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP1074 10040
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
 Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
 Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
 Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

TidyBlot WESTERN BLOT DETECTION REAGENT:HRP (STAR209P)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M401395:220715'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)