

Datasheet: AHP1168

BATCH NUMBER 160802

Description:	RABBIT ANTI CD5L (C-TERMINAL)
Specificity:	CD5L (C-TERMINAL)
Other names:	SP-ALPHA
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			0.5 - 1.0ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species	Human
Species Cross Reactivity	<p>Reacts with: Mouse</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
Product Form	Purified IgG - liquid

Antiserum Preparation Antisera to human CD5L were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.02% Sodium Azide

Approx. Protein Concentrations IgG concentration 1.0mg/ml

Immunogen 13 amino acid peptide sequence near the carboxy terminus of human CD5L.

External Database Links

UniProt:
[O43866](#) [Related reagents](#)

Entrez Gene:
[922](#) CD5L [Related reagents](#)

Synonyms API6

RRID AB_808449

Specificity **Rabbit anti Human CD5L antibody** detects the C-terminal (CT) region of human CD5 like (CD5L), also known as Apoptosis Inhibitor expressed by Macrophages (AIM), and SP alpha. CD5L is a soluble ~38 kDa glycoprotein produced by mature macrophages, and has been shown to support their survival and enhance phagocytic function. CD5L expression is upregulated in response to LXR/RXR heterodimer activation.

CD5L has been shown to mediate the survival of macrophages in atherosclerotic plaques.

Western Blotting AHP1168 detects a band of approximately 38kDa in Raji cell lysates.

References 1. Li, Y. *et al.* (2008) Extracellular Nampt promotes macrophage survival via a nonenzymatic interleukin-6/STAT3 signaling mechanism. [J Biol Chem. 283 \(50\): 34833-43.](#)

Further Reading 1. Sarrias, M.R. *et al.* (2005) A role for human Sp alpha as a pattern recognition receptor. [J Biol Chem. 280 \(42\): 35391-8.](#)
2. Arai, S. *et al.* (2005) A role for the apoptosis inhibitory factor AIM/Spalpha/Api6 in atherosclerosis development. [Cell Metab. 1 \(3\): 201-13.](#)

Storage Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at:
<https://www.bio-rad-antibodies.com/SDS/AHP1168>
10040

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)
[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M363827:200529'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)