

Datasheet: AHP1184

BATCH NUMBER 163824

Description:	RABBIT ANTI SUMO-1 (N-TERMINAL)
Specificity:	SUMO-1 (N-TERMINAL)
Other names:	SMALL UBIQUITIN-RELATED MODIFIER-1
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			0.5 - 1.0ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species

Human

Species Cross Reactivity

Reacts with: Rat, Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Antiserum Preparation

Antisera to human SUMO-1 were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.02% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	14 amino acid peptide sequence near the amino terminus of human SUMO-1.
External Database Links	<p>UniProt: P63165 Related reagents</p> <p>Entrez Gene: 7341 SUMO1 Related reagents</p>
Synonyms	SMT3C, SMT3H3, UBL1
RRID	AB_2198094
Specificity	<p>Rabbit anti Human SUMO-1 antibody recognizes the N-terminal (NT) region of small ubiquitin-related modifier-1 (SUMO-1). SUMO-1 is an ~11 kDa protein which structurally related to ubiquitin, and is added to proteins as a form of post-translational modification known as sumoylation. Sumoylation can be cell cycle regulated, stimulated upon stress or induced by DNA damage, and binding elicits individual effects for each protein. Effects include changes to protein stability, nuclear-cytosolic transport, and transcriptional regulation (mostly transcriptional repression).</p> <p>Alternate transcriptional splice variants encoding different isoforms have been characterized.</p>
Western Blotting	AHP1184 detects a band of approximately 15kDa in HL-60 cell lysates.
Further Reading	<ol style="list-style-type: none"> Melchior, F. & Hengst, L. (2002) SUMO-1 and p53. Cell Cycle. 1 (4): 245-9. Verger, A. <i>et al.</i> (2003) Modification with SUMO. A role in transcriptional regulation. EMBO Rep. 4 (2): 137-42.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP1184

10040

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)

[FITC](#)

Sheep Anti Rabbit IgG (STAR35...)

[RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...)

[Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR36...)

[DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

Goat Anti Rabbit IgG (H/L) (STAR124...)

[HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M381647:210512'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)