

Datasheet: AHP1189

Description:	RABBIT ANTI BMI-1
Specificity:	BMI-1
Other names:	POLYCOMB GROUP RING FINGER PROTEIN 4
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			0.5 - 1.0ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
Species Cross Reactivity	Reacts with: Rat, Mouse N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to human BMI-1 were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.02% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	15 amino acid peptide sequence near the centre of human BMI-1.

**External Database
Links**

UniProt:

[P35226](#) [Related reagents](#)

Entrez Gene:

[648](#) BMI1 [Related reagents](#)

Synonyms

PCGF4, RNF51

RRID

AB_2065518

Specificity

Rabbit anti BMI-1 antibody recognizes 37 kDa polycomb complex protein BMI-1, involved in segment specification, cell growth and maintenance, transcriptional regulation and chromatin modification. BMI-1 is a proto-oncogene which cooperates with the myc oncogene to produce B-lymphomas. It is required in the regulation of normal adult and leukaemic stem cells, making BMI-1 a potential target for stem cell proliferation and anti-cancer therapies. BMI-1 represses transcription at the ink4 locus which contains tumour suppressors p16 and p19Arf.

BMI-1 is a component of the Polycomb group (PcG) multiprotein PRC1 complex which ubiquitinates histone-2A at lys119, resulting in a heritable change in the ability of chromatin to express genes. BMI-1 is also required to stimulate the E3 ubiquitin ligase activity of RNF2/RING2 within the PRC1 complex.

The Notch ligand DLL4 (delta-like protein 4) enhances T cell proliferation, by directly influencing BMI-1 transcription.

Further Reading

1. Raaphorst, F. *et al.* (2003) Self-renewal of hematopoietic and leukemic stem cells: a central role for the Polycomb-group gene Bmi-1. [Trends Immunol. 10:522-4.](#)
2. Alkema, M.J. *et al.* (1993) Characterization and chromosomal localization of the human proto-oncogene BMI-1. [Hum Mol Genet. 2 \(10\): 1597-603.](#)
3. Lessard, J. & Sauvageau, G. (2003) Bmi-1 determines the proliferative capacity of normal and leukaemic stem cells. [Nature. 423 \(6937\): 255-60.](#)

Storage

Store at +4°C or at -20°C if preferred.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10040 available at:
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M363847:200529'

Printed on 11 Aug 2020

© 2020 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)