

## Datasheet: AHP1212

**BATCH NUMBER 147097**

<b>Description:</b>	GOAT ANTI HUMAN TNF ALPHA
<b>Specificity:</b>	TNF ALPHA
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			2.5 - 5.0ug/ml
ELISA	▪			0.25 - 1.0ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1 - 0.2ug/ml
Functional Assays	▪			0.012 - 0.014ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

**(1) This product heat-mediated antigen retrieval prior to staining paraffin sections. Citrate buffer, pH6.0, is recommended for this purpose.**

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - lyophilised
<b>Reconstitution</b>	Reconstitute with 0.1ml distilled water. For long term storage the addition of 0.09% sodium azide is recommended. NB: For functional studies do not add azide.
<b>Antiserum Preparation</b>	Antisera to human TNF $\alpha$ were raised by repeated immunisations of goats with highly purified antigen. Purified IgG prepared by affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0mg/ml after reconstitution
Immunogen	<a href="#">Recombinant human TNF<math>\alpha</math></a>
External Database Links	<p><b>UniProt:</b>  <a href="#">P01375</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">7124</a>    TNF    <a href="#">Related reagents</a></p>
Synonyms	TNFA, TNFSF2
RRID	AB_872048
Specificity	<p><b>Goat anti Human TNF alpha antibody</b> recognizes human Tumour Necrosis Factor alpha (TNF-<math>\alpha</math>), a 17.4kDa multi-functional pro-inflammatory cytokine primarily secreted by macrophages, but also produced by T and B lymphocytes, activated monocytes and fibroblasts.</p> <p>TNF-<math>\alpha</math> exists in both a transmembrane and mature soluble form and has many biological properties, including activation of transcription factor NF-<math>\kappa</math>B, the induction of apoptosis in various tumour cell lines and the stimulation of cell differentiation and proliferation under certain conditions. TNF-<math>\alpha</math> also plays a role in inflammation, septic shock, autoimmune disease and rheumatoid arthritis.</p>
ELISA	This purified human TNF $\alpha$ antibody may be used in an indirect ELISA or as the capture reagent in a sandwich ELISA with our <a href="#">biotinylated TNF<math>\alpha</math> antibody</a> (AHP1212B) as the detection reagent and <a href="#">recombinant human TNF<math>\alpha</math></a> (PHP051) as the standard
Histology Positive Control Tissue	Human Breast Carcinoma
References	1. van Horssen, R. <i>et al.</i> (2006) TNF- $\alpha$ in cancer treatment: Molecular insights, antitumor effects, and clinical utility. <a href="#">Oncologist. 11: 397-408.</a>
Storage	<p>Prior to reconstitution store at -20°C.  After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody.</p>

Should this product contain a precipitate we recommend microcentrifugation before use.

---

<b>Guarantee</b>	Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.
------------------	---

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10162 available at: <a href="https://www.bio-rad-antibodies.com/SDS/AHP1212">https://www.bio-rad-antibodies.com/SDS/AHP1212</a> 10162
--------------------------------------	--

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

### Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M372473:200706'

Printed on 23 May 2025

---

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)