

Datasheet: AHP1248

BATCH NUMBER 153784

Description:	RABBIT ANTI HUMAN TAP1 (C-TERMINAL)
Specificity:	TAP1 (C-TERMINAL)
Other names:	TRANSPORTER ASSOCIATED WITH ANTIGEN PROCESSING 1
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/10000
Immunoprecipitation	▪			
Western Blotting	▪			1/1000
Functional Assays			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to TAP1 (C-terminal) were raised by repeated immunisations of rabbits with highly purified antigen.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Sodium Azide (NaN ₃)

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	TAP1 peptide corresponding to the C-terminal of the human protein conjugated to Keyhole Limpet Hemocyanin (KLH).
External Database Links	<p>UniProt: Q03518 Related reagents</p> <p>Entrez Gene: 6890 TAP1 Related reagents</p>
Synonyms	ABCB2, PSF1, RING4, Y3
RRID	AB_2303049
Specificity	<p>Rabbit anti Human TAP1 antibody detects an epitope within the C-terminal region of human TAP1, also known as transporter associated with antigen processing 1. TAP1 is an ~81 kDa membrane-associated protein which is a member of the ATP-binding cassette (ABC) transporter superfamily and together with the homologous subunit TAP2 forms the functional TAP complex.</p> <p>TAP1 is involved in the transport of antigenic peptides from the cytosol into the endoplasmic reticulum, where they associate with major histocompatibility complex (MHC) class I molecules. These complexes are presented at the cell surface to T-lymphocytes, which identify and remove infected or malignantly transformed cells.</p> <p>Rabbit anti Human TAP1 antibody (AHP1248) does not cross react with human TAP2.</p>
Western Blotting	AHP1248 detects a band of approximately 72kDa
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP1248 10040</p>
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M363869:200529'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)