

Datasheet: AHP1278

**BATCH NUMBER 164604**

<b>Description:</b>	RABBIT ANTI HUMAN CTGF
<b>Specificity:</b>	CTGF
<b>Other names:</b>	CONNECTIVE TISSUE GROWTH FACTOR
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			
ELISA	▪			0.5 - 2.0ug/ml
Western Blotting	▪			0.1 - 0.2ug/ml
Immunofluorescence	▪			
Functional Assays	▪			6.0 - 9.0ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

**(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.**

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - lyophilized
<b>Reconstitution</b>	Reconstitute with 0.1ml distilled water. For long term storage the addition of 0.09% sodium azide is recommended. NB: For functional studies do not add azide.

**Antiserum Preparation** Antisera to human CTGF were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	None present
<b>Carrier Free</b>	Yes
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0mg/ml after reconstitution.
<b>Immunogen</b>	<a href="#">Recombinant human CTGF</a> .
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P29279</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">1490</a>    CTGF    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	CCN2, HCS24, IGFBP8
<b>RRID</b>	AB_2087370
<b>Specificity</b>	<p><b>Rabbit anti Human CTGF antibody</b> recognizes human CTGF (connective tissue growth factor), a 323 amino acid ~38 kDa polypeptide member of the CCN family of secreted cysteine rich regulatory proteins. CTGF is produced by vascular endothelial and umbilical vein cells, and interacts with several other growth factors, including transforming growth factor beta (TGF-beta), vascular endothelial growth factor (VEGF) and bone morphogenetic proteins (BMPs).</p> <p>CTGF is a heparin-binding pleiotropic growth factor involved in chondrogenesis and potentially skeletogenesis. CTGF promotes the adhesion of epithelial cells and fibroblasts, stimulates integrin expression, and also promotes endothelial cell survival, migration and adhesion, implicating CTGF in angiogenesis. CTGF induces extracellular matrix formation, including the deposition of excessive collagen in fibrotic diseases, and has also been implicated in wound healing.</p>
<b>ELISA</b>	This purified human CTGF antibody may ne used in an indirect ELISA or as the capture reagent in a sandwich ELISA with our <a href="#">biotinylated human CTGF antibody</a> (AHP1278B) as the detection reagent and <a href="#">recombinant human CTGF</a> (PHP174) as the standard.
<b>References</b>	<ol style="list-style-type: none"> <li>Westermann, D. <i>et al.</i> (2009) Gene deletion of the kinin receptor B1 attenuates cardiac inflammation and fibrosis during the development of experimental diabetic cardiomyopathy. <a href="#">Diabetes. 58: 1373-81.</a></li> <li>Waddell, J.M. <i>et al.</i> (2011) CTGF expression is up-regulated by PROK1 in early pregnancy and influences HTR-8/Svneo cell adhesion and network formation. <a href="#">Hum Reprod. 26 (1): 67-75.</a></li> <li>Vasilieva. O.V. <i>et al.</i> (2016) Connective tissue growth factor (CTGF) in the human dermis through ontogenesis <a href="#">Russian Journal of Developmental Biology. 47 (2): 63-68.</a></li> </ol>

4. Zarobkiewicz, M.K. *et al.* (2018) Decrease in Lipid Droplets in Adrenal Cortex of Male Wistar Rats after Chronic Exposure to Energy Drinks. [Medicina \(Kaunas\). 54 \(5\):90.](#)

---

**Further Reading** 1. Moussad, E.E. & Brigstock, D.R. (2000) Connective tissue growth factor: what's in a name? [Mol Genet Metab. 71 \(1-2\): 276-92.](#)

---

**Storage** Prior to reconstitution store at -20°C.  
After reconstitution store at -20°C.

This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee** 12 months from date of despatch

---

**Health And Safety Information** Material Safety Datasheet documentation #10294 available at: <https://www.bio-rad-antibodies.com/SDS/AHP1278>  
10294

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)  
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)  
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)  
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

### Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)  
[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M399197:220628'

Printed on 18 Jan 2024

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)