

Datasheet: AHP1282

BATCH NUMBER 155895

Description:	RABBIT ANTI HUMAN IGFBP-7
Specificity:	IGFBP-7
Other names:	INSULIN-LIKE GROWTH FACTOR BINDING PROTEIN-7
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.5 - 2.0ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1 - 0.2ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - lyophilised
Reconstitution	Reconstitute with 0.1ml distilled water. For long term storage the addition of 0.09% sodium azide is recommended.
Antiserum Preparation	Antisera to human IGFBP-7 were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative	None present

Stabilisers

Carrier Free Yes

Approx. Protein Concentrations IgG concentration 1.0mg/ml after reconstitution.

Immunogen [Recombinant human IGFBP-7](#)

External Database Links**UniProt:**[Q16270](#) [Related reagents](#)**Entrez Gene:**[3490](#) IGFBP7 [Related reagents](#)

Synonyms MAC25, PSF

RRID AB_2248672

Specificity **Rabbit anti Human IGFBP-7 antibody** recognizes human IGFBP-7 (insulin-like growth factor binding protein-7), also known as mac25, a 256 amino acid secreted protein and member of a growing family of IGFBPs which regulate the biological activities and distribution of the insulin-like growth factors (IGFs), through binding to IGF-I and IGF-II.

IGFBP-7 is widely expressed throughout normal tissues, but is also expressed in many tumour types at an increased or reduced level, and acts as a suppressor of breast and prostate cancer cell lines. IGFBP-7 inhibits IGF-induced differentiation of skeletal muscle (myogenesis) and has implications in bone formation, by its possible effects on the actions of parathyroid hormone (PTH) and cortisol.

ELISA This purified human IGFBP-7 antibody may be used in an ELISA [recombinant human IGFBP-7](#) (PHP178) as the standard.

Further Reading

1. Ruan, W.J. *et al.* (2006) IGFBP7 plays a potential tumor suppressor role against colorectal carcinogenesis with its expression associated with DNA hypomethylation of exon 1. [J Zhejiang Univ Sci B. 7 \(11\): 929-32.](#)
2. Pereira, R.C. & Canalis, E. (1999) Parathyroid hormone increases mac25/insulin-like growth factor-binding protein-related protein-1 expression in cultured osteoblasts. [Endocrinology. 140 \(5\): 1998-2003.](#)

Storage Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.

This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 6 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10162 available at:
<https://www.bio-rad-antibodies.com/SDS/AHP1282>
10162

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

 Email: antibody_sales_us@bio-rad.com

Worldwide

 Tel: +44 (0)1865 852 700

 Fax: +44 (0)1865 852 739

 Email: antibody_sales_uk@bio-rad.com

Europe

 Tel: +49 (0) 89 8090 95 21

 Fax: +49 (0) 89 8090 95 50

 Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M373919:201013'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)