

Datasheet: AHP1283

BATCH NUMBER 165794

Description:	RABBIT ANTI HUMAN INTERLEUKIN-17F
Specificity:	IL-17F
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.5 - 2.0ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1 - 0.2ug/ml
Functional Assays	▪			0.9 - 1.2 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - lyophilized
Reconstitution	Reconstitute with 0.1ml distilled water. For long term storage the addition of 0.09% sodium azide is recommended. N.B. For functional studies do not add sodium azide.
Antiserum Preparation	Antisera to human IL-17F were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG was prepared by affinity chromatography.
Buffer Solution	Phosphate buffered saline

Preservative Stabilisers	None present
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	Recombinant human IL-17F.
External Database Links	<p>UniProt: Q96PD4 Related reagents</p> <p>Entrez Gene: 112744 IL17F Related reagents</p>
Synonyms	IL24
RRID	AB_2125420
Specificity	<p>Rabbit anti Human interleukin-17F antibody recognizes human IL-17F (interleukin-17F), a disulphide-linked homodimer consisting of two 133 amino acid chains, which belongs to the IL-17 gene family and shares the highest sequence homology with IL-17A.</p> <p>IL-17F is expressed by activated CD4+ T-cells and monocytes, and acts as a stimulator of cytokine production, including granulocyte colony-stimulating factor (GM-CSF), IL-6 and IL-8. IL-7F also stimulates the production of IL-2, MCP-1 (monocyte chemotactic protein-1) and TGFB (Transforming growth factor beta) by endothelial cells, and inhibits endothelial cell angiogenesis.</p> <p>Rabbit anti Human interleukin-17F antibody (AHP1283) has been used in neutralisation assays where 50% inhibition (ND50) of the biological activity of 50ng/ml of human IL-17F was achieved using 0.9-1.2µg/ml of antibody.</p>
References	1. Starnes, T. <i>et al.</i> (2001) Cutting edge: IL-17F, a novel cytokine selectively expressed in activated T cells and monocytes, regulates angiogenesis and endothelial cell cytokine production. J Immunol. 167 (8): 4137-40.
Storage	<p>Prior to reconstitution store at -20°C. After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10294 available at:
<https://www.bio-rad-antibodies.com/SDS/AHP1283>
10294

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M399236:220628'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)