

Datasheet: AHP1440

BATCH NUMBER 157986

Description:	RABBIT ANTI HUMAN PTK7 (N-TERMINAL)
Specificity:	PTK7 (N-TERMINAL)
Other names:	CCK-4
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			0.5 - 1.0ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species

Human

Species Cross Reactivity

Reacts with: Rat, Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation	Antisera to human PTK7 (NT) were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.02% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	A 15 amino acid peptide located near the N-terminal region of human PTK7.
External Database Links	<p>UniProt: Q13308 Related reagents</p> <p>Entrez Gene: 5754 PTK7 Related reagents</p>
Synonyms	CCK4
RRID	AB_2174361
Specificity	<p>Rabbit anti Human PTK7 antibody recognizes the N-terminal region of human Protein tyrosine kinase-like 7 (PTK7), a member of the receptor tyrosine kinase family, otherwise known as Colon Carcinoma Kinase 4 (CCK4). PTK7 is highly expressed in the kidney, liver, lung, pancreas and melanocytes with weaker expression being detected in the thyroid gland, ovaries, brain, heart and skeletal muscle.</p> <p>PTK7 may play a role in embryonic development. It is thought to function as a cell adhesion molecule and as a novel regulator of the planar cell polarity (PCP) pathway in vertebrates.</p> <p>PTK7 may be involved in the tumour progression. A loss of PTK7 expression has been observed in melanoma cell lines, whereas increased levels of PTK7 have been detected in high-throughput acute myeloid leukaemia samples. In addition, PTK7 is not expressed in normal colon but has been detected in colon carcinomas and as such is a potential marker for tumour progression.</p>
Western Blotting	AHP1440 detects a band of approximately 86kDa in colon tissue lysates.
Further Reading	<ol style="list-style-type: none"> Lu, X. <i>et al.</i> (2004) PTK7/CCK-4 is a novel regulator of planar cell polarity in vertebrates. Nature. 430 (6995): 93-8. Easty, D.J. <i>et al.</i> (1997) Loss of expression of receptor tyrosine kinase family genes PTK7 and SEK in metastatic melanoma. Int J Cancer. 71 (6): 1061-5.
Storage	Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP1440 10040
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Recommended Useful Reagents

- [ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)
[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M363930:200529'

Printed on 18 Jan 2024