

Datasheet: AHP1624

**BATCH NUMBER 168779**

<b>Description:</b>	RABBIT ANTI TEM7
<b>Specificity:</b>	TEM7
<b>Other names:</b>	TEM3
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin (1)	▪			2.5ug/ml
Western Blotting	▪			0.5 - 1.0ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

**(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.**

### Target Species

Human

### Species Cross Reactivity

Reacts with: Mouse, Rat

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG - liquid

### Antiserum Preparation

Antiserum to human TEM7 was raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.02% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0mg/ml
<b>Immunogen</b>	A peptide corresponding to a 16 amino acid sequence from near the center of human TEM7.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">Q8IUK5</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">57125</a>    PLXDC1    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	TEM3, TEM7
<b>RRID</b>	AB_2166226
<b>Specificity</b>	<p><b>Rabbit anti Human TEM7 antibody</b> recognises human Tumour endothelial marker 7, also known as TEM7, a 55.8 kDa protein belonging to the plexin family.</p> <p>TEM7 is strongly expressed in endothelial cells from colorectal cancer and from primary cancers of the lung, liver, brain, breast and pancreas. High expression of TEM7 is associated with metastasis and poor survival of patients with osteogenic sarcoma.</p> <p>There appear to be at least four isoforms of TEM7 which include membrane bound, secreted and intracellular forms. TEM7 is not detectable in endothelial cells from normal tissues.</p>
<b>Histology Positive Control Tissue</b>	Mouse liver
<b>Western Blotting</b>	AHP1624 detects a band of approximately 56kDa in rat liver cell lysates.
<b>Storage</b>	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/SDS/AHP1624">https://www.bio-rad-antibodies.com/SDS/AHP1624</a>

---

**Regulatory**For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)

[FITC](#)

Sheep Anti Rabbit IgG (STAR35...)

[RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...)

[Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR36...)

[DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

Goat Anti Rabbit IgG (H/L) (STAR124...)

[HRP](#)

### Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)**North & South** Tel: +1 800 265 7376**Worldwide**

Tel: +44 (0)1865 852 700

**Europe**

Tel: +49 (0) 89 8090 95 21

**America** Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M382815:210513'

**Printed on 29 Aug 2024**