

Datasheet: AHP1757

**BATCH NUMBER 163944**

<b>Description:</b>	GOAT ANTI HUMAN TRANSFERRIN
<b>Specificity:</b>	TRANSFERRIN
<b>Format:</b>	Serum
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Turbidimetry	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
-----------------------	-------

<b>Product Form</b>	Serum - liquid
---------------------	----------------

<b>Antiserum Preparation</b>	Antisera to Human Transferrin were raised by repeated immunisation of goats with highly purified antigen.
------------------------------	---

<b>Buffer Solution</b>	TRIS buffered saline
------------------------	----------------------

<b>Preservative Stabilisers</b>	0.1% Sodium Azide (NaN <sub>3</sub> )
---------------------------------	---------------------------------------

<b>Immunogen</b>	Purified human Transferrin
------------------	----------------------------

### External Database Links

#### UniProt:

[P02787](#)

[Related reagents](#)

#### Entrez Gene:

---

<b>RRID</b>	AB_2202010
-------------	------------

---

<b>Specificity</b>	<p><b>Goat anti human transferrin polyclonal antibody</b> recognizes human Transferrin, also known as Serotransferrin, an iron binding glycoprotein expressed by the liver and secreted into plasma.</p> <p>Transferrin has an important role in transporting iron from the intestine, reticuloendothelial system, and liver parenchymal cells to all proliferating cells in the body.</p> <p>This reagent is available in bulk quantities and may be made available in other formats. Please enquire for further details.</p>
--------------------	--

---

<b>Storage</b>	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
----------------	---

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10342 available at: <a href="https://www.bio-rad-antibodies.com/SDS/AHP1757">https://www.bio-rad-antibodies.com/SDS/AHP1757</a> 10342
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M382584:210513'

Printed on 15 Mar 2024