

Datasheet: AHP1757

BATCH NUMBER 168934

Description:	GOAT ANTI HUMAN TRANSFERRIN
Specificity:	TRANSFERRIN
Format:	Serum
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Turbidimetry	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Serum - liquid
Antiserum Preparation	Antisera to Human Transferrin were raised by repeated immunisation of goats with highly purified antigen.
Buffer Solution	TRIS buffered saline
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Immunogen	Purified human Transferrin

External Database Links

UniProt:

[P02787](http://www.uniprot.org/entry/P02787)

[Related reagents](#)

Entrez Gene:

RRID	AB_2202010
Specificity	Goat anti human transferrin polyclonal antibody recognizes human Transferrin, also known as Serotransferrin, an iron binding glycoprotein expressed by the liver and secreted into plasma. Transferrin has an important role in transporting iron from the intestine, reticuloendothelial system, and liver parenchymal cells to all proliferating cells in the body. This reagent is available in bulk quantities and may be made available in other formats. Please enquire for further details.
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10342 available at: https://www.bio-rad-antibodies.com/SDS/AHP175710342
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M429055:240315'

Printed on 29 Aug 2024