

Datasheet: AHP1757

Description:	GOAT ANTI HUMAN TRANSFERRIN
Specificity:	TRANSFERRIN
Format:	Serum
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Turbidimetry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
-----------------------	-------

Product Form	Serum - liquid
---------------------	----------------

Antiserum Preparation	Antisera to Human Transferrin were raised by repeated immunisation of goats with highly purified antigen.
------------------------------	---

Buffer Solution	TRIS buffered saline
------------------------	----------------------

Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
---------------------------------	---------------------------------------

Immunogen	Purified human Transferrin
------------------	----------------------------

External Database Links

UniProt:

[P02787](#) [Related reagents](#)

Entrez Gene:

[7018](#) TF [Related reagents](#)

RRID AB_2202010

Specificity **Goat anti human transferrin polyclonal antibody** recognizes human Transferrin, also known as Serotransferrin, an iron binding glycoprotein expressed by the liver and secreted into plasma.

Transferrin has an important role in transporting iron from the intestine, reticuloendothelial system, and liver parenchymal cells to all proliferating cells in the body.

This reagent is available in bulk quantities and may be made available in other formats. Please enquire for further details.

Storage Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10342 available at:
10342: <https://www.bio-rad-antibodies.com/uploads/MSDS/10342.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M364022:200529'

Printed on 22 Mar 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)