

Datasheet: AHP1765

BATCH NUMBER 160921

Description:	GOAT ANTI HUMAN KAPPA LIGHT CHAIN
Specificity:	KAPPA LIGHT CHAIN (BOUND AND FREE)
Format:	Serum
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Turbidimetry	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
-----------------------	-------

Product Form	Serum - liquid
---------------------	----------------

Antiserum Preparation	Antisera to Human Kappa Light chains were raised by repeated immunisation of goats with highly purified antigen.
------------------------------	------------------------------------------------------------------------------------------------------------------

Buffer Solution	TRIS buffered saline
------------------------	----------------------

Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
---------------------------------	----------------------------------------

Immunogen	Purified human kappa light chains
------------------	-----------------------------------

External Database Links

UniProt:
[P01834](#) [Related reagents](#)

Entrez Gene:

RRID	AB_1605166
-------------	------------

Specificity	Goat anti Human kappa light chain antibody recognises the Human kappa immunoglobulin light chain which comprises approximately 60% of light chains in the human.
--------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10057 available at: https://www.bio-rad-antibodies.com/SDS/AHP176510057
--------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M364028:200529'

Printed on 12 Aug 2023