

Datasheet: AHP1765

Description:	GOAT ANTI HUMAN KAPPA LIGHT CHAIN
Specificity:	KAPPA LIGHT CHAIN (BOUND AND FREE)
Format:	Serum
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Turbidimetry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
-----------------------	-------

Product Form	Serum - liquid
---------------------	----------------

Antiserum Preparation Antisera to Human Kappa Light chains were raised by repeated immunisation of goats with highly purified antigen.

Buffer Solution	TRIS buffered saline
------------------------	----------------------

Preservative Stabilisers	<0.1% Sodium Azide (NaN ₃)
---------------------------------	--

Immunogen	Purified human kappa light chains
------------------	-----------------------------------

External Database Links

UniProt:

[P01834](#) [Related reagents](#)

Entrez Gene:

[3514](#) IGKC [Related reagents](#)

RRID	AB_1605166
Specificity	Goat anti Human kappa light chain antibody recognises the Human kappa immunoglobulin light chain which comprises approximately 60% of light chains in the human.
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10057 available at: 10057: https://www.bio-rad-antibodies.com/uploads/MSDS/10057.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M389164:210806'

Printed on 29 Aug 2021