

Datasheet: AHP1809

**BATCH NUMBER 163811**

<b>Description:</b>	GOAT ANTI HUMAN CD194
<b>Specificity:</b>	CD194
<b>Other names:</b>	CCR4
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			1/50 - 1/100
ELISA	▪			1/10,000 - 1/30,000
Immunoprecipitation			▪	
Western Blotting	▪			1/100 - 1/500

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Antiserum Preparation</b>	Antiserum to human CD194 was raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared by affinity chromatography.
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative</b>	0.1% Sodium Azide (NaN <sub>3</sub> )
<b>Stabilisers</b>	0.1% Bovine Serum Albumin
<b>Approx. Protein</b>	IgG concentration 1.0mg/ml

## Concentrations

---

**Immunogen** Synthetic peptide SNYYLYESIPKPCTKEGIKAFGE corresponding to amino acids 17-39 of human CD194.

---

## External Database

### Links

**UniProt:**[P51679](#)[Related reagents](#)**Entrez Gene:**[1233](#)

CCR4

[Related reagents](#)

---

**Synonyms**

CMKBR4

---

**RRID**

AB\_2074495

---

**Specificity**

**Goat anti Human CD194 antibody** recognises human CD194, otherwise known as CCR4 or C-C chemokine receptor type 4. CD194 is a transmembrane glycoprotein and member of the G-protein coupled receptor 1 family, preferentially expressed by CD4+ Th2 (type 2 helper T) cells, which acts as a high affinity receptor for the chemokines MDC (macrophage-derived chemokine), MCP-1, MIP-1, RANTES and TARC.

---

**Flow Cytometry**

Use 10ul of the suggested working dilution to label  $1 \times 10^6$  cells in 100ul.

---

**Histology Positive Control Tissue**

Spleen

---

**Storage**

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety Information**

Material Safety Datasheet documentation #10041 available at:  
10041: <https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf>

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

### Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

**North & South** Tel: +1 800 265 7376

**Worldwide** Tel: +44 (0)1865 852 700

**Europe** Tel: +49 (0) 89 8090 95 21

**America** Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To  
find a

batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M382957:210513'

**Printed on 08 Feb 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)