

Datasheet: AHP1831

BATCH NUMBER 163901

Description:	RABBIT ANTI HUMAN BDNF
Specificity:	BDNF
Other names:	BRAIN-DERIVED NEUROTROPHIC FACTOR
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.5 - 2.0ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1 - 0.2ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - lyophilized
Reconstitution	<p>Reconstitute with 0.1ml distilled water</p> <p>Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended.</p>

Antiserum Preparation Antiserum to human BDNF were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0mg/ml after reconstitution
Immunogen	Recombinant human BDNF .
External Database Links	<p>UniProt: P23560 Related reagents</p> <p>Entrez Gene: 627 BDNF Related reagents</p>
RRID	AB_2064219
Specificity	<p>Rabbit anti Human BDNF antibody recognizes human brain-derived neurotrophic factor (BDNF), also known as abrineurin. BDNF is a secreted 119 amino acid ~30 kDa protein with an additional 18 amino acid signal peptide and a 110 amino acid propeptide region. BDNF is a potent neurotrophin and member of the NGF-beta family, secreted from dense core vesicles, which promotes the growth, differentiation and survival of neurons within, or directly connected to, the central nervous system (CNS).</p> <p>BDNF signals through binding to the high affinity receptor trkB, and also binds to the low affinity neurotrophin receptor LNGFR. BDNF is involved in the modulation of short-term synaptic function and in the formation of long-term memory (Alonso et al. 2002). BDNF is modulated during progression of Alzheimer's with an early increase and reduced levels in patients with latter stages of Alzheimer's (Laske et al. 2005). Decreased levels of BDNF are also noted in the progression of Huntington's disease (Zuccato et al. 2008) .</p>
ELISA	This purified human BDNF antibody may be used in an indirect ELISA or as the capture reagent in a sandwich ELISA with our biotinylated human BDNF antibody (AHP1831B) as the detection reagent and recombinant human BDNF (PHP213) as the standard.
Storage	<p>Prior to reconstitution store at -20°C. After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety	Material Safety Datasheet documentation #10294 available at:

Information <https://www.bio-rad-antibodies.com/SDS/AHP183110294>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M402202:220718'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)