

Datasheet: AHP1833

BATCH NUMBER 149319

Description:	RABBIT ANTI HUMAN CD120b
Specificity:	CD120b
Other names:	TNF-R2
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.5 - 2.0ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1 - 0.2ug/ml
Immunofluorescence	▪			3.0ug/ml
Functional Assays	▪			15.0 - 20.0ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - lyophilised
Reconstitution	<p>Reconstitute with 0.1ml distilled water</p> <p>Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended.</p> <p>N.B. For functional studies do not add sodium azide</p>

Antiserum Preparation Antiserum to human CD120b were raised by repeated immunisation of rabbits with highly

purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative Stabilisers	None present
---------------------------------	--------------

Carrier Free	Yes
---------------------	-----

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml after reconstitution
---------------------------------------	--------------------------------------------------

Immunogen	Recombinant human CD120b (PHP215)
------------------	---------------------------------------------------

External Database Links	UniProt: P20333 Related reagents Entrez Gene: 7133 TNFRSF1B Related reagents
--------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Synonyms	TNFBR, TNFR2
-----------------	--------------

RRID	AB_2240832
-------------	------------

Specificity	<p>Rabbit anti human CD120b recognizes the soluble form of human CD120b, otherwise known as sTNF receptor type II (sTNFR1I), a member of the tumour necrosis factor receptor (TNFR) superfamily, produced by proteolytic cleavage of the membrane form by the metalloprotease TACE.</p> <p>CD120b plays a key role in regulating the inflammatory response, acting as a receptor for both tumour necrosis factor alpha (TNF-alpha) and TNF-beta, and is also involved in programmed cell death and anti-viral activity. sTNFR1I acts as an inhibitor of TNF-alpha induced responses, by acting as a decoy receptor.</p>
--------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

ELISA	This purified human CD120b antibody may be used in an indirect ELISA or as the capture reagent in a sandwich ELISA with a biotinylated human CD120b antibody (AHP1833B) as the detection reagent and recombinant human CD120b (PHP215) as the standard.
--------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Western Blotting	This purified human CD120b antibody may be used in Western Blotting applications under either reducing or non-reducing conditions with recombinant human CD120b (PHP215) as the positive control.
-------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Storage	Prior to reconstitution store at -20°C. After reconstitution store at -20°C.
----------------	---------------------------------------------------------------------------------

This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee	6 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10162 available at: https://www.bio-rad-antibodies.com/SDS/AHP1833 10162
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M372950:200730'

Printed on 04 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)