

Datasheet: AHP1870

Description:	SHEEP ANTI GUINEA PIG CALRETININ
Specificity:	CALRETININ
Other names:	CALB2
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			1/500
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			1/1000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Guinea Pig
-----------------------	------------

Species Cross Reactivity	<p>Reacts with: Rat</p> <p>Based on sequence similarity, is expected to react with: Human, Mouse, Primate</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
---------------------------------	---

Product Form	Purified IgG - liquid
---------------------	-----------------------

Antiserum Preparation	Antisera to Calretinin were raised by repeated immunisations of sheep with highly purified antigen. Purified IgG was prepared from whole serum by Protein A chromatography.
------------------------------	---

Buffer Solution	TRIS buffered saline
Preservative Stabilisers	None present
Immunogen	Purified guinea pig calretinin
External Database Links	<p>UniProt:</p> <p>P47728 Related reagents</p> <p>P22676 Related reagents</p> <p>Q08331 Related reagents</p> <p>Entrez Gene:</p> <p>117059 Calb2 Related reagents</p> <p>794 CALB2 Related reagents</p> <p>12308 Calb2 Related reagents</p>
Synonyms	CAB29
RRID	AB_1720089
Specificity	<p>Sheep anti calretinin antibody recognizes calretinin, a ~29 kDa calcium binding protein found mostly in the central nervous system and retina. It is a cytoplasmic protein involved in the regulation of calcium concentration, which is important in the regulation of a large variety of neuronal functions, including neurotransmitter release, ionic channel permeability, enzyme activity, and gene transcription. Mice deficient in calretinin show impairment in motor coordination.</p> <p>Calretinin is a sensitive and specific marker for mast cells.</p>
Western Blotting	AHP1870 detects a band of approximately 29kDa in rat cerebellum.
Storage	<p>Store at -20°C only.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #10170 available at: https://www.bio-rad-antibodies.com/SDS/AHP1870</p> <p>10170</p>
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Sheep IgG (H/L) (5184-2304...) [Biotin](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M428168:240301'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)