

Datasheet: AHP2174

Description:	RABBIT ANTI HUMAN CD294 (N-TERMINAL)
Specificity:	CD294 (N-TERMINAL)
Other names:	Putative G-protein coupled receptor 44
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

RRID AB_10709827

Applications This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin	▪			5ug/ml
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting	▪			
Functional Assays			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Human

Species Cross Reactivity Based on sequence similarity, is expected to react with: Mouse, Rat
N.B. Antibody reactivity and working conditions may vary between species.

Product Form Purified IgG - liquid

Antiserum Preparation Antisera to human CD294 were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG prepared from whole serum by affinity chromatography.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.02% Sodium Azide (NaN₃)

Approx. Protein Concentrations IgG concentration 1.0mg/ml

Immunogen	A synthetic peptide corresponding to an 18 amino acid sequence from near the amino terminus of human CD294 – KLH conjugated.
External Database Links	<p>UniProt: Q9Y5Y4 Related reagents Q9Z2J6 Related reagents</p> <p>Entrez Gene: 11251 GPR44 Related reagents 14764 Gpr44 Related reagents</p>
Synonyms	Crth2, CRTH2, DL1R
Specificity	<p>Rabbit anti Human CD294 antibody recognizes CD294, also known as chemoattractant receptor homologous molecule expressed on Th2 cells (CRTH2) and GRP44. CD294 is a G protein-coupled receptor which has a high affinity for prostaglandin D2 (PGD2).</p> <p>CD294 is predominantly expressed on Th2 cells, basophils and eosinophils, and its activation by PDG2 mediates chemotaxis of these cells to inflammatory sites. CD294 may be a useful marker for Th2 cells, and enhances their function and signalling with a possible role in Th2 cell differentiation.</p> <p>PGD2 is thought to contribute to bronchoconstriction, eosinophilia and mucus production in allergic asthma, and CD294 may contribute to this disease.</p>
Immunohistology	<p>*This product does not require protein digestion pre-treatment of paraffin sections e.g. trypsin or pronase. This product requires antigen retrieval using heat treatment prior to staining of paraffin sections.</p> <p>Sodium citrate buffer pH 6.0 is recommended for this purpose.</p>
Histology Positive Control Tissue	Human adrenal gland
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Sheep Anti Rabbit IgG (STAR54...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Recommended Useful Reagents

[HISTAR DETECTION SYSTEM \(STAR3000A\)](#)

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025C\)](#)

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M343808:190111'

Printed on 20 May 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)