

Datasheet: AHP2188

BATCH NUMBER 166076

Description:	RABBIT ANTI HUMAN WNT-3a
Specificity:	WNT-3a
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			
ELISA	▪			0.5 - 2.0ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1 - 0.2ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species	Human
Product Form	Purified IgG - lyophilized
Reconstitution	Reconstitute with 0.1ml distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended.

Antiserum Preparation Antiserum to human WNT-3a were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0mg/ml after reconstitution
Immunogen	Recombinant human WNT-3a.
External Database Links	<p>UniProt: P56704 Related reagents</p> <p>Entrez Gene: 89780 WNT3A Related reagents</p>
RRID	AB_10844208
Specificity	<p>Rabbit anti Human WNT-3a antibody recognizes human WNT-3a, a secreted protein and member of the Wnt family of signalling proteins, acting as a ligand for members of the frizzled family of seven transmembrane receptors. Wnt-3a plays key roles in many developmental and cell fate decisions, including embryonic patterning, CNS development, and cytoskeletal formation.</p> <p>Wnt-3a acts as a promoter of the proliferation of ventral midbrain dopaminergic neuron precursor cells, and has been implicated in the progressive differentiation of anterior neural cells, as well as regulating the survival and expansion of neural progenitors derived from human embryonic stem cells. Wnt3a acts upstream of the Notch signalling pathway, in the regulation of somitogenesis, controlled by a segmentation clock, consisting of molecular oscillators in the Notch, Wnt3a and Fgf8 pathways(Aulehla <i>et al.</i> 2003).</p>
ELISA	AHP2188 may be used in an indirect ELISA or as a capture reagent in a sandwich ELISA with (AHP2188B) as the detection antibody.
Further Reading	<ol style="list-style-type: none"> Lai, Y.H. <i>et al.</i> (2011) SOX4 interacts with plakoglobin in a Wnt3a-dependent manner in prostate cancer cells. BMC Cell Biol. 12: 50. Yan, L. <i>et al.</i> (2011) Activation of the canonical Wnt/β-catenin pathway in ATF3-induced mammary tumors. PLoS One. 6:e16515.
Storage	<p>Prior to reconstitution store at -20°C. After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10294 available at:
<https://www.bio-rad-antibodies.com/SDS/AHP2188>
10294

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M402166:220718'

Printed on 19 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)