

Datasheet: AHP2224

BATCH NUMBER 180201

Description:	RABBIT ANTI VPS34
Specificity:	VPS34
Other names:	PI3-KINASE TYPE 3
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			2.0 - 4.0ug/ml
Immunofluorescence	▪			1/100

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Human

Species Cross Reactivity Reacts with: Mouse, Rat
N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form Purified IgG - liquid

Antiserum Preparation Antiserum to human VSP34 was raised by repeated immunisation of rabbits with highly

purified antigen. Purified IgG was prepared by affinity chromatography.

Buffer Solution	Phosphate buffered saline.
------------------------	----------------------------

Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
---------------------------------	--

Approx. Protein Concentrations	IgG concentration 1.0mg/ml
---------------------------------------	----------------------------

Immunogen	KLH conjugated synthetic peptide corresponding to amino acids 22-34 of human VPS34, The corresponding sequence is identical in mouse and rat.
------------------	---

External Database Links	UniProt: Q8NEB9 Related reagents O88763 Related reagents Q6PF93 Related reagents Entrez Gene: 5289 PIK3C3 Related reagents 65052 Pik3c3 Related reagents 225326 Pik3c3 Related reagents
--------------------------------	--

Synonyms	Vps34, VPS34
-----------------	--------------

RRID	AB_10845944
-------------	-------------

Specificity	<p>Rabbit anti VPS34 antibody recognises Phosphatidylinositol 3-kinase catalytic subunit type 3, also known as VPS34 (vacuolar protein sorting 34), PI3-kinase type 3 or PIK3C3. VPS34 is a ubiquitously expressed class III PI3K (phosphoinositide 3-kinase), and member of the PI3/PI4-kinase family, which utilizes phosphatidylinositol (PtdIn) as a substrate, and functions as the catalytic subunit of the PI3K complex.</p> <p>VPS34 plays a role in the transport of lysosomal enzyme precursors to lysosomes, and in endocytosis and vesicular trafficking, and has been strongly implicated as a regulator of autophagy, forming a multi-protein complex with the proautophagic tumor suppressors Beclin1/Atg6, Bif-1, and UVRAG, in the initiation of autophagosome formation. Distinct complexes of VPS34 regulates endocytic processes vital for late-stage autophagosome-lysosome fusion.</p>
--------------------	---

Western Blotting	AHP224 detects a band of approximately 105kDa in HEK-293T cell lysates.
-------------------------	---

Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	---

Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP2224 10040
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M364153:200529'

Printed on 25 Mar 2023