

Datasheet: AHP2259GA

BATCH NUMBER 149979

Description:	RABBIT ANTI BOVINE NEUROFILAMENT H 200kDa
Specificity:	NEUROFILAMENT H 200kDa
Other names:	Neurofilament heavy polypeptide
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			1.0ug/ml
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			0.1ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species	Bovine
Species Cross Reactivity	<p>Reacts with: Rat, Human, Mouse</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
Product Form	Purified IgG - liquid

Antiserum Preparation Antiserum to bovine NFH was raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared by affinity chromatography.

Buffer Solution TRIS-glycine buffered saline, NaCl

Preservative Stabilisers 0.05% Sodium Azide (NaN₃)

Approx. Protein Concentrations IgG concentration 1.0mg/ml

Immunogen GST-tagged bovine NFH recombinant protein.

External Database Links

UniProt:

[P12036](#) [Related reagents](#)

[P19246](#) [Related reagents](#)

[P16884](#) [Related reagents](#)

Entrez Gene:

[4744](#) NEFH [Related reagents](#)

[380684](#) Nefh [Related reagents](#)

[24587](#) Nefh [Related reagents](#)

Synonyms Kiaa0845, KIAA0845, Nfh, NFH

Specificity **Rabbit anti Bovine Neurofilament H 200 kDa antibody** recognizes bovine intermediate 200 kDa Neurofilament H (NFH), a 1026 amino acid ~200 kDa protein involved in the maintenance of neuronal integrity.

Neurofilaments are a major component of the cellular cytoskeleton, acting as the most abundant support for the axon cytoplasm. They exist as three neurofilament proteins: 68/70 kDa (NFL) 160 kDa (NFM) and 200 kDa (NFH), which are expressed in both the central and peripheral nervous systems. Although typically expressed by neurons, neurofilaments are also expressed by neuroblastomas, neuromas, gangliogliomas, ganglioneuromas, and ganglioneuroblastomas, as well as paragangliomas, carcinoids, skin neuroendocrine carcinomas, and oat cell carcinomas of the lung ([Lehto et al. 1983](#)).

Defects in NFH are responsible for susceptibility to the neurodegenerative disorder amyotrophic lateral sclerosis (ALS) which affects upper and lower motor neurons, and results in fatal paralysis ([Mendonça et al. 2005](#)).

Histology Positive Control Tissue Bovine or rat cerebellum

Western Blotting AHP2259GA detects a band of approximately 200kDa in bovine cerebellum tissue lysates.

References 1. Leypoldt, F. *et al.* (2002) Neuronal differentiation of cultured human NTERA-2cl.D1 cells leads to increased expression of synapsins. [Neurosci Lett. 324 \(1\): 37-40.](#)

2. Ignatov, A. *et al.* (2003) Role of the G-protein-coupled receptor GPR12 as high-affinity receptor for sphingosylphosphorylcholine and its expression and function in brain development. [J Neurosci. 23: 907-14.](#)
3. Kaneko, S. *et al.* (2006) Protecting axonal degeneration by increasing nicotinamide adenine dinucleotide levels in experimental autoimmune encephalomyelitis models. [J Neurosci. 26: 9794-804.](#)
4. Li, Y. *et al.* (2008) Transplanted olfactory ensheathing cells incorporated into the optic nerve head ensheath the retinal ganglion cell axons: possible relevance to glaucoma. [Neurosci Lett. 440: 251-4.](#)
5. Zhang, Z. *et al.* (2011) The role of single cell derived vascular resident endothelial progenitor cells in the enhancement of vascularization in scaffold-based skin regeneration. [Biomaterials. 32: 4109-17.](#)
6. Chen, K. *et al.* (2016) Differential Histopathological and Behavioral Outcomes Eight Weeks after Rat Spinal Cord Injury by Contusion, Dislocation, and Distraction Mechanisms. [J Neurotrauma. 33 \(18\): 1667-84.](#)

Further Reading 1. Liu Q *et al.* (2004) Neurofilament proteins in neurodegenerative diseases. [Cell Mol Life Sci. 61 \(24\): 3057-75.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10511 available at: <https://www.bio-rad-antibodies.com/SDS/AHP2259GA>
10511

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M387455:210629'

Printed on 21 Mar 2025

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)