

Datasheet: AHP2268

BATCH NUMBER 163962

Description:	GOAT ANTI WNT-3
Specificity:	WNT-3
Other names:	INT4
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			3.0 - 5.0ug/ml
ELISA	▪			1:16000
Immunoprecipitation			▪	
Western Blotting	▪			1.0 - 3.0ug/ml
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1)This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species	Human
Species Cross Reactivity	<p>Reacts with: Mouse</p> <p>Based on sequence similarity, is expected to react with: Dog, Rat</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>

Product Form	Purified IgG - liquid
Antiserum Preparation	Antiserum to Wnt-3 was raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared by affinity chromatography using the immunizing peptide
Buffer Solution	TRIS buffered saline
Preservative	0.02% Sodium Azide (NaN ₃)
Stabilisers	0.5% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Synthetic peptide sequence C-GRGHNTRTEKRKEK from the internal region of Wnt-3 (NP_110380.1).
External Database Links	<p>UniProt: P56703 Related reagents</p> <p>Entrez Gene: 7473 WNT3 Related reagents</p>
Synonyms	INT4
Specificity	<p>Goat anti Wnt-3 antibody recognises both human and mouse proto-oncogene Wnt-3, a member of the Wnt family, and ligand of the frizzled family of seven transmembrane receptors.</p> <p>Wnt-3, otherwise known as Int-4, is involved in cell fate decisions and patterning during embryogenesis, cell signalling during neural tube development, and implicated in the tumorigenesis of several cancers including breast, lung, hepatocellular(HCC), and gastric carcinoma.</p> <p>Over-expression of Wnt-3 plays a pivotal role in the resistance of HER2-overexpressing breast tumors to the breast cancer treatment drug Trastuzumab, also known as Herceptin (Wu et al. 2012). The effects of Wnt-3 over-expression included the transactivated expression of EGFR (epidermal growth factor receptor) an increased expression of B-catenin within the nucleus, an increase in N-cadherin and the neural crest transcription factor SLUG, as well as promoting an EMT (epithelial to mesenchymal transition) response.</p> <p>Goat anti Wnt-3 polyclonal antibody exhibits strong fibrous staining in myoepithelial cells of human breast tissue in paraffin sections.</p> <p>Goat anti Wnt-3 polyclonal antibody stains mouse hippocampal primary astrocytes in ICC/IF (Okamoto et al. 2011).</p>

Histology Positive Control Tissue	Human breast
Western Blotting	AHP2268 detects a band of approximately 40kDa in HEK293 cell lysates
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10058 available at: https://www.bio-rad-antibodies.com/SDS/AHP2268 10058
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M382847:210513'

Printed on 01 Mar 2024